KOMBITEST™ CD3 FITC / CD19 PE

Cat. No. / Kat. č.: ED7009

ENGLISH

1. Specification

Content	50 tests, 1 ml	
Usage	20 µl per test	
Specificity	CD3	CD19
Clone	MEM-57	LT19
Isotype (mouse)	IgG2a	IgG1
Fluorochrome	FITC	PE
λ excitation	488 nm	488 nm
Emission maximum	525 nm	575 nm

2. Intended use

The KOMBITEST™ CD3 FITC / CD19 PE is designed for percentage enumeration of mature human T and B lymphocytes in erythrocyte-lysed whole blood.

3. Principle

This test is based on specific binding of monoclonal antibodies to the antigenic determinants expressed on the surface of leukocytes. The monoclonal antibodies are labeled with different fluorochromes which are excited via laser beam from a flow cytometer during analysis. Subsequent emissions of light from fluorochromes of each cell are collected and analyzed by a flow cytometer. The fluorescence intensity differences enable separation of cell subsets based on expression of analyzed antigens.

Specific staining of blood cells is performed by incubation of blood samples with the reagent followed by a lysis of red blood cells. Afterwards, unaffected leukocytes are subjected to analysis by a flow cytometer.

4. Specificity

The monoclonal antibody MEM-57 reacts with gamma-epsilon and delta-epsilon dimers of human CD3 complex, a part of a bigger multisubunit T cell receptor complex (CD3/TCR) expressed on peripheral blood T lymphocytes and mature thymocytes. The monoclonal antibody MEM-57 was assigned to CD3 during the Human Leukocyte Differentiation Antigen workshop (HLDA4 WS Code: T 96).

The antibody LT19 reacts with CD19 (B4), a 95 kDa type I transmembrane glycoprotein of immunoglobulin superfamily,

expressed on B lymphocytes and follicular dendritic cells; lost on plasma cells.

5. Reagent provided

The reagent contains premixed combination of mouse monoclonal antibody against human CD3 antigen (clone MEM-57) labeled with Fluorescein isothiocyanate (FITC), and mouse monoclonal antibody against human CD19 antigen (clone LT19) labeled with R-phycoerythrin (PE). The labeled antibodies are diluted at optimum concentration in phosphate buffered saline (PBS) containing 15mM sodium azide and 0.2% (w/v) high-grade protease free Bovine Serum Albumin (BSA) as a stabilizing agent. The content of a vial (1 ml) is sufficient for 50 tests.

6. Storage

Store vial in the dark at 2-8°C. Do not freeze.

7. Precautions

- Intended for professional use only.
- Do not use after expiration date stamped on vial label.
- Avoid prolonged exposure to light.
- The content of the vial must not freeze.
- Avoid contamination of the reagent.
- Any non-performance of staining protocol may produce false results.
- The reagent contains sodium azide (NaN₃) which is highly toxic in pure form. However, the concentration in the reagent (15mM) is not considered as hazardous. When disposing the reagent, flush the sink with large volume of water to avoid accumulation of explosive metal-azide in plumbing.
- Blood samples are considered as potentially infectious and must be handled with care. Avoid all contact of the sample with the skin, eyes and mucosa.

8. Necessary material not supplied

Material necessary for collection of peripheral blood, test tubes for staining of blood samples (e.g. 12 × 75 mm), automatic pipettes with disposable tips, vortex mixer, centrifuge, commercial lysing solution, phosphate buffered saline (PBS), FITC/PE Negative Control (Cat. No. ED7032), flow cytometer.

9. Staining protocol

1. Collect peripheral blood in a sterile tube with an anticoagulant (e.g. Heparin, EDTA).
2. Add 20 µl of KOMBITEST™ CD3 FITC / CD19 PE reagent to a test tube, and the same amount of FITC/PE Negative Control to a control tube.
3. Add 100 µl of blood sample to each tube. Vortex the tubes.
4. Incubate tubes for 20-30 minutes at room temperature in the dark.
5. Perform lysis of red cells using lysing solution. It is recommended to use a commercial lysing solution containing formaldehyde as a fixative. Follow the instructions of the lysing solution manufacturer.
6. Centrifuge tubes for 5 minutes at 300 g.
7. Remove supernatant and resuspend pellet with 3-4 ml of PBS.
8. Centrifuge tubes for 5 minutes at 300 g.
9. Remove supernatant and resuspend pellet with 0.3 – 0.5 ml of PBS.
10. Analyze samples immediately using flow cytometer or store samples at 2-8°C in the dark and analyze within 24 hours provided that cells were fixed.

10. Data analysis

Analyze sample stained with KOMBITEST™ CD3 FITC / CD19 PE using a flow cytometer. Make compensation of fluorescent signals prior or after data acquisition. Visualize compensated data on the side-scatter (SSC) versus forward-scatter (FSC) plot. Set the gate for lymphocyte population as shown in figure 1. Alternatively, set the optimal lymphocyte gate using KOMBITEST™ CD45 FITC / CD14 PE (refer the datasheet for lymphocyte gate assessment procedure). Then make a CD3 FITC versus CD19 PE dot-plot of lymphocyte population as shown in figure 2. Separate populations using appropriate gate and calculate the percentage of T lymphocytes situated in lower-right quadrant (CD3+CD19-) and B lymphocytes situated in upper-left quadrant (CD3-CD19+) on the dot-plot. The regions corresponding to the negative populations should be set up using control cells which were stained by FITC/PE Negative Control.

11. Expected values

Results obtained in different laboratories may vary. Each laboratory should establish a normal range of cell subsets using its own test conditions. In our laboratory, the KOMBITEST™ CD3 FITC / CD19 PE was tested on 40 blood samples of healthy people. Obtained results are given in the table below.

Parameter	Mean (%)	SD	CV (%)
B lymphocytes	8.4	4.3	51.2
T lymphocytes	72.2	7.7	10.7

12. Limitations

- Flow cytometer may produce false results if the device has not been aligned and maintained appropriately.

- Data may be incorrectly interpreted if fluorescent signals were compensated wrongly or if gates were positioned inaccurately.
- Blood samples from abnormal patients may exhibit abnormal values of positive cells.
- Red blood cells from abnormal patients may be resistant to lysis using lysing solutions.
- In case of hyperleukocytose sample, it is recommended to dilute blood sample with PBS to obtain leukocyte density approximately 5×10^6 leukocytes/ml.
- Blood samples should be stained and analyzed within 24 hours from the blood collection.

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
Fig. 1: Delimitation of lymphocyte population

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Alarcon B et al. (2006) T-cell antigen-receptor stoichiometry: pre-clustering for sensitivity. *EMBO Rep.* 7: 490-5

Kuhns MS et al. (2006) Deconstructing the form and function of the TCR/CD3 complex. *Immunity* 24:133-9

van Zelm MC et al. (2006) An antibody-deficiency syndrome due to mutations in the CD19 gene. *N Engl J Med.* 354: 1901-1912

Lin CW et al. (2005) CD94 1A transcripts characterize lymphoblastic lymphoma/leukemia of immature natural killer cell origin with distinct clinical features. *Blood.* 106: 3567-74

Huang Y and Wange RL (2004) T cell receptor signaling: beyond complex complexes. *J Biol Chem.* 279: 28827-28830

Elias F et al. (2003) Strong cytosine-guanosine-independent immunostimulation in humans and other primates by synthetic oligodeoxynucleotides with PyNTTTTGT motifs. *J Immunol.* 171: 3697-704

Inabe K and Kurosaki T (2002) Tyrosine phosphorylation of B-cell adaptor for phosphoinositide 3-kinase is required for Akt activation in response to CD19 engagement. *Blood* 99: 584-589

Fujimoto M et al. (1999) CD19 amplifies B lymphocyte signal transduction by regulating Src-Family protein tyrosine kinase activation. *J Immunol.* 162: 7088-7094


Hilgert I. et al. (1993) Therapeutic in vivo use of the A1-CD3 monoclonal antibody. *Transplantation* 55: 435

Horejsi V. et al. (1988) Monoclonal antibodies against human leucocyte antigens. II. Antibodies against CD45 (T200), CD3 (T3), CD43, CD10 (CALLA), transferrin receptor (T9), a novel broadly expressed 18-kDa antigen (MEM-43) and a novel antigen of restricted expression (MEM-74). *Folia Biol. (Praha)* 34: 23

Leukocyte Typing IV., Knapp W. et al. (Eds.), Oxford University Press (1989); p. 293.

Leukocyte Typing III., McMichael M.J. et al. (Eds.), Oxford University Press (1987); p.611.

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