EXCELLYSE I

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ENGLISH

1. Intended use

The EXCELLYSE I lysing solution permits red blood cell lysis following antibody staining of human peripheral blood leukocytes. This solution is optimized for use with BD Biosciences flow cytometers. Other manufacturer's flow cytometer may also be used without considerable sample quality loss.

EXCELLYSE I lysing solution is appropriate for use with EXBIO single colour monoclonal antibodies and KOMBITEST™ reagents and may be used in both lyse/wash and lyse/no wash protocol.

2. Introduction

Leukocyte analysis and detection in peripheral blood requires elimination of interfering cells, mainly erythrocytes. So far, Ficoll density gradient method was used to separate leukocytes from whole blood. This method is rather time consuming and may lead to a loss of certain lymphocyte subsets. Direct blood sample staining followed by red blood cell lysis therefore takes place in clinical laboratories as a fast and accurate method for whole blood flow cytometry analysis. EXCELLYSE I lysing solution provides for excellent results. Erythrocytes are lysed with only minimum debris left and leukocytes are fixed keeping cell subsets well distinguished and keeping cell markers preserved.

3. Reagent provided

EXCELLYSE I lysing solution is a ready to use reagent provided as 30 ml or 100 ml package of proprietary buffered solution containing < 5 % formaldehyde. This quantity is sufficient for 150 or 500 tests respectively.

4. Storage

Store the EXCELLYSE I lysing solution at room temperature.

5. Precautions

- For research use only.
- Do not use after expiration date stamped on vial label.
- Do not use if any discoloration or precipitation occurs.
- Do not freeze.
- Avoid contamination of the reagent.
- Any non-performance of lysing protocol may produce false results.

- The reagent contains < 5% formaldehyde which is toxic and a suspected carcinogen. It may be allergenic. Avoid any contact with skin, eyes, mucous membranes and clothing. Appropriate care should be taken.
- Blood samples are considered as potentially infectious and must be handled with care. Avoid all contact of the sample with the skin, eyes and mucosa.

6. Necessary material not supplied

Material necessary for collection of peripheral blood, test tubes for blood sample staining (e.g. 12 × 75 mm), automatic pipettes with disposable tips, vortex mixer, centrifuge, phosphate buffered saline (PBS), deionized water, appropriate fluorescent-dye-labeled primary/secondary antibody, flow cytometer.

7. Lyse/wash blood lysing protocol

1. Collect peripheral blood in a sterile tube with an anticoagulant (e.g. Heparin, EDTA).
2. Follow manufacturer's instructions for whole blood antibody staining.
3. Add 200 µl of EXCELLYSE I lysing solution per 100 µl of whole blood. Mix the tube with a vortex mixer.
4. Incubate the tube for 10 minutes at room temperature in the dark.
5. Add 3-4 ml of deionized water to the tube and incubate for about 5-10 minutes, until the blurry blood sample solution gets clear.
6. Centrifuge the tube for 5 minutes at 300 g.
7. Remove supernatant and resuspend pellet with 0.3 – 0.5 ml of PBS.
8. Analyze the sample immediately using flow cytometer or store the sample at 2-8 °C in the dark and analyze within 24 hours. No further cell fixation is required.

8. Lyse/no wash blood lysing protocol

1. Collect peripheral blood in a sterile tube with an anticoagulant (e.g. Heparin, EDTA).
2. Follow instructions for whole blood antibody staining.
3. Add 200 µl of EXCELLYSE I lysing solution per 100 µl of whole blood. Mix the tube with a vortex mixer.
4. Incubate the tube for 10 minutes at room temperature in the dark.
5. Add 1 ml of deionized water to the tube, mix well, and incubate for about 5-10 minutes, until the blurry blood sample solution gets clear.
6. Analyze the sample immediately using flow cytometer or store the sample at 2-8°C in the dark and analyze within 24 hours. No further cell fixation is required.

9. Limitations

- Flow cytometer may produce false results if the device has not been aligned and maintained appropriately.
- Data may be incorrectly interpreted if fluorescent signals were compensated wrongly or if gates were positioned inaccurately.
- Blood samples from abnormal patients may exhibit abnormal values of positive cells.
- Red blood cells from abnormal patients may be resistant to lysis using lysing solutions.
- In case of hyperleukocytose sample, it is recommended to dilute blood sample with PBS to obtain leukocyte density approximately 5×10^6 leukocytes/ml.
- Blood samples should be stained and analyzed within 24 hours from the blood collection.

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
Fig. 1: Peripheral blood leukocytes dot-plot from lysed whole blood using EXCELLYSE I lysing solution, analyzed on BD FACSCanto™ cytometer.

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
Fig. 2: Peripheral blood leukocytes dot-plot from lysed whole blood using EXCELLYSE I lysing solution, analyzed on Dako Cyan™ cytometer.

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
Fig. 3: Peripheral blood leukocytes dot-plot from lysed whole blood using EXCELLYSE I lysing solution, analyzed on Beckman Coulter FC-500 cytometer.

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Trademarks / Registrované názvy

BD FACSCanto™ is a registered trademark of BD Biosciences and Company.

Dako Cyan™ is a registered trademark of Beckman Coulter Company.