

# Monoclonal Antibody to CD14, PE conjugated (CD14 PE)

Cat.No. ED7015

## 1. Intended purpose

The reagent CD14 PE permits identification and enumeration of cell populations expressing human CD14 antigen in whole blood using flow cytometry.

## 2. Test principle

This test is based on specific binding of monoclonal antibody to the antigenic determinant expressed on the surface of leukocytes. The monoclonal antibody is labeled with fluorochrome which is excited via laser beam from a flow cytometer during analysis. Subsequent emission of light from fluorochromes of each cell is collected and analyzed by a flow cytometer. The fluorescence intensity differences enable separation of cell subsets based on expression of analyzed antigen. Specific staining of blood cells is performed by incubation of blood samples with the reagent followed by a lysis of red blood cells. Afterwards, unaffected leukocytes are subjected to analysis by a flow cytometer.

## 3. Reagents provided

The reagent contains mouse monoclonal antibody against human CD14 antigen (clone MEM-15) which was purified by affinity chromatography and labeled with R-Phycoerythrin (PE). The labeled antibody is diluted at optimum concentration in stabilizing phosphate buffered saline (PBS) solution containing 15mM sodium azide. The content of a vial (2 ml) is sufficient for 100 tests.

## Product specification

|                  |                 |
|------------------|-----------------|
| Content          | 100 tests, 2 ml |
| Usage            | 20 µl per test  |
| Specificity      | Human CD14      |
| Clone            | MEM-15          |
| Isotype          | Mouse IgG1      |
| Fluorochrome     | PE              |
| λ excitation     | 488 nm          |
| Emission maximum | 575 nm          |

## 4. Materials required but not provided

Test tubes for blood staining (e.g. 12 × 75 mm)  
Commercial lysing solution  
Phosphate buffered saline (PBS)  
Isotype control antibody (mouse IgG1 PE)

## 5. Equipment required

Automatic pipettes with disposable tips  
Vortex mixer  
Centrifuge  
Flow cytometer with excitation laser 488 nm and proper filters

## 6. Storage and handling

Store the vial at 2 - 8 °C. Keep away from sunlight. Do not freeze. Do not aliquot. Expiration date is stated on a vial label and on outer packaging.

## 7. Warnings, precautions and limitations of use

- Intended for In Vitro Diagnostic use in laboratories outside USA and Canada. This reagent is in conformity with the European Directive 98/79/EC.
- Do not use reagent after expiration date.
- Avoid reagents contamination.
- Avoid prolonged exposure to light.
- The content of the vial must not freeze.
- Any non-performance of staining protocol may produce false results.
- The reagent contains sodium azide (Na<sub>2</sub>N<sub>3</sub>) which is highly toxic in pure form. However, the concentration in the reagent (15mM) is not considered as hazardous. When disposing the reagent, flush the sink with a large volume of water.
- Blood samples are considered as potentially infectious and must be handled with care. Avoid all contact of the sample with the skin, eyes and mucosa.
- In case of hyperleukocytose sample, it is recommended to dilute blood sample with PBS to obtain leukocyte density approximately  $5 \times 10^5$  leukocytes/ml.
- Blood samples from abnormal patients may exhibit abnormal values of positive cells.
- Data may be incorrectly interpreted if fluorescent signals were compensated

wrongly or if gates were positioned inaccurately.

- Flow cytometer may produce false results if the device has not been aligned and maintained appropriately.
- Red blood cells from abnormal patients may be resistant to lysis using lysing solutions.
- Blood samples should be stained and analyzed within 24 hours from the blood collection.

## 8. Specimen

Use the peripheral human blood in a sterile tube with an anticoagulant (Heparin or EDTA). Blood must be stored at room temperature.

## 9. Procedure


### Staining protocol

- Add 20 µl of CD14 PE reagent to a test tube, and the necessary amount of isotype control to a control tube.
- Add 100 µl of blood sample to each tube. Vortex the tubes.
- Incubate tubes for 20 - 30 minutes at room temperature in the dark.
- Perform lysis of red cells using lysing solution. It is recommended to use a commercial lysing solution containing formaldehyde as a fixative. Follow the instructions of the lysing solution manufacturer.
- Centrifuge tubes for 5 minutes at 300 g.
- Remove supernatant and resuspend pellet with 3 - 4 ml of PBS.
- Centrifuge tubes for 5 minutes at 300 g.
- Remove supernatant and resuspend pellet with 0.3 - 0.5 ml of PBS.
- Analyze samples immediately using flow cytometer or store samples at 2 - 8 °C in the dark and analyze within 24 hours provided that cells were fixed.

### Flow Cytometric Analysis

Analyze sample stained with CD14 PE using a flow cytometer. Visualize recorded data using appropriate plot such as side-scatter (SSC) versus PE intensity as shown in figure 1.

Fig. 1: Leukocytes stained with CD14 PE reagent



The brightest population (CD14+) belongs to monocytes. Set suitable gates for analysis. The region corresponding to the negative population should be set up using control cells which were stained by isotype control antibody.

## 10. Analytical performance

### Specificity

The antibody MEM-15 reacts with CD14, a 53-55 kDa GPI (glycosylphosphatidylinositol) linked membrane glycoprotein expressed on monocytes, macrophages and weakly on granulocytes; it is expressed by most tissue macrophages. This antibody also reacts with soluble forms of CD14 found in serum and in the urine of some nephrotic patients. The monoclonal antibody MEM-15 was assigned to CD14 during the Human Leukocyte Differentiation Antigen workshop (HLDA3 WS Code: M 252).

### 11. Clinical performance

#### Expected values

Results obtained in different laboratories may vary. Each laboratory should establish a normal range of cell subsets using its own test conditions.

### 12. References

- Asai Y et al. (2007) Soluble CD14 Discriminates Slight Structural Differences between Lipid As That Lead to Distinct Host Cell Activation. *J Immunol.* 179: 7674-83
- Lodrup Carlsen KC and Granum B (2007) Soluble CD14: role in atopic disease and recurrent infections, including otitis media. *Curr Allergy Asthma Rep.* 7: 436-43
- Fernández-Real JM et al. (2003) CD14 monocyte receptor, involved in the inflammatory cascade, and insulin sensitivity. *J Clin Endocrinol Metab.* 88: 1780-4

Juan TS et al. (1995) Identification of a domain in soluble CD14 essential for lipopolysaccharide (LPS) signaling but not LPS binding. *J Biol Chem.* 270: 17237-42

Bazil V et al. (1986) Biochemical characterization of a soluble form of the 53-kDa monocyte surface antigen. *Eur J Immunol.* 16: 1583-9

Leukocyte Typing VI., Kishimoto T. et al. (Eds.), Garland Publishing Inc. (1997).

Leukocyte Typing V., Schlossman S. et al. (Eds.), Oxford University Press (1995).

Leukocyte Typing IV., Knapp W. et al. (Eds.), Oxford University Press (1989).

Leukocyte Typing III., McMichael A. J. et al (Eds.), Oxford University Press (1987).

## 13. Manufacturer

EXBIO Praha, a.s.  
Nad Safinou II 341  
25250 Vestec  
Czech Republic

[info@exbio.cz](mailto:info@exbio.cz)  
[technical@exbio.cz](mailto:technical@exbio.cz)  
[orders@exbio.cz](mailto:orders@exbio.cz)  
[www.exbio.cz](http://www.exbio.cz)

## 14. Trademarks

N/A

## 15. Revision History

- Version 1, ED7015\_IFU\_v1  
Initial Release
- Version 2, ED7015\_IFU\_v2  
Merging three language mutations into one document. The address was changed: "Nad Safinou II 341".
- Version 3, ED7015\_IFU\_v3  
Precautions section was changed - "Intended for professional use only." - removed. "Intended for In Vitro Diagnostic use in laboratories outside USA and Canada. This CE-IVD reagent is in conformity with the European In Vitro Diagnostic Medical Device Directive 98/79/EC." - added.
- Version 4, ED7015\_IFU\_v4  
Reagent provided section was changed: text "stabilizing" added, "solution" - added and "0.2% (w/v) high-grade protease free Bovine Serum Albumin (BSA) as a stabilizing agent" - removed.
- Version 5, ED7015\_IFU\_v5  
The company logo changed. IFU layout changed. "Keep away from sunlight." - added. Postal code changed: "25250 Vestec"

# exbio

---

## Monoclonal Antibody to CD14, PE conjugated (CD14 PE)

100 tests | Cat.No. ED7015



### Instructions for Use

Version: ED7015\_IFU\_v5\_EN  
Date of Issue: 29-07-2020

EN

### Symbols

|   |                                    |
|---|------------------------------------|
| The REF symbol consists of the letters 'REF' inside a rectangular border.     | Catalogue number                   |
| The LOT symbol consists of the letters 'LOT' inside a rectangular border.     | Batch code                         |
| The use-by date symbol is an hourglass icon.                                  | Use-by date                        |
| The temperature limits symbol is a thermometer icon.                          | Temperature limits                 |
| The keep away from sunlight symbol is a sun with a slash through it.          | Keep away from sunlight            |
| The IVD symbol consists of the letters 'IVD' inside a rectangular border.     | In vitro diagnostic medical device |
| The CE marking symbol consists of the letters 'C' and 'E' in a stylized font. | CE marking of conformity           |
| The consult instructions for use symbol is an icon of an open book.           | Consult instructions for use       |
| The manufacturer symbol is an icon of a factory.                              | Manufacturer                       |

The product is intended for In Vitro Diagnostic Use. In vivo diagnostic or therapeutic applications are strictly forbidden.

Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at [www.exbio.cz](http://www.exbio.cz). EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s..