

Monoclonal Antibody to CD56, PE conjugated (CD56 PE)

Cat.No. ED7019

1. Intended purpose

The reagent CD56 PE permits identification and enumeration of cell populations expressing human CD56 antigen in whole blood using flow cytometry.

2. Test principle

This test is based on specific binding of monoclonal antibody to the antigenic determinant expressed on the surface of leukocytes. The monoclonal antibody is labeled with fluorochrome which is excited via laser beam from a flow cytometer during analysis. Subsequent emission of light from fluorochromes of each cell is collected and analyzed by a flow cytometer. The fluorescence intensity differences enable separation of cell subsets based on expression of analyzed antigen. Specific staining of blood cells is performed by incubation of blood samples with the reagent followed by a lysis of red blood cells. Afterwards, unaffected leukocytes are subjected to analysis by a flow cytometer.

3. Reagents provided

The reagent contains mouse monoclonal antibody against human CD56 antigen (clone MEM-188) which was purified by affinity chromatography and labeled with R-Phycoerythrin (PE). The labeled antibody is diluted at optimum concentration in stabilizing phosphate buffered saline (PBS) solution containing 15mM sodium azide. The content of a vial (2 ml) is sufficient for 100 tests.

Product specification

Content	100 tests, 2 ml
Usage	20 µl per test
Specificity	Human CD56
Clone	MEM-188
Isotype	Mouse IgG2a
Fluorochrome	PE
λ excitation	488 nm
Emission maximum	575 nm

4. Materials required but not provided

Test tubes for blood staining (e.g. 12 × 75 mm)
Commercial lysing solution
Phosphate buffered saline (PBS)
Isotype control antibody (mouse IgG2a PE)

5. Equipment required

Automatic pipettes with disposable tips
Vortex mixer
Centrifuge
Flow cytometer with excitation laser 488 nm and proper filters

6. Storage and handling

Store the vial at 2-8°C. Keep away from sunlight. Do not freeze. Do not aliquot. Expiration date is stated on a vial label and on outer packaging.

7. Warnings, precautions and limitations of use

- Intended for In Vitro Diagnostic use in laboratories outside USA and Canada. This reagent is in conformity with the European Directive 98/79/EC.
- Do not use reagent after expiration date.
- Avoid reagent contamination.
- Avoid prolonged exposure to light.
- The content of the vial must not freeze.
- Any non-performance of staining protocol may produce false results.
- The reagent contains sodium azide (NaN₃) which is highly toxic in pure form. However, the concentration in the reagent (15mM) is not considered as hazardous. When disposing the reagent, flush the sink with a large volume of water.
- Blood samples are considered as potentially infectious and must be handled with care. Avoid all contact of the sample with the skin, eyes and mucosa.
- In case of hyperleukocytose sample, it is recommended to dilute blood sample with PBS to obtain leukocyte density approximately 5 × 10⁶ leukocytes/ml.
- Blood samples from abnormal patients may exhibit abnormal values of positive cells.
- Data may be incorrectly interpreted if fluorescent signals were compensated wrongly or if gates were positioned

inaccurately.

- Flow cytometer may produce false results if the device has not been aligned and maintained appropriately.
- Red blood cells from abnormal patients may be resistant to lysis using lysing solutions.
- Blood samples should be stained and analyzed within 24 hours from the blood collection.

8. Specimen

Use the peripheral human blood in a sterile tube with an anticoagulant (Heparin or EDTA). Blood must be stored at room temperature.

9. Procedure

Staining protocol

- Add 20 µl of CD56 PE reagent to a test tube, and the necessary amount of isotype control to a control tube.
- Add 100 µl of blood sample to each tube. Vortex the tubes.
- Incubate tubes for 20 - 30 minutes at room temperature in the dark.
- Perform lysis of red cells using lysing solution. It is recommended to use a commercial lysing solution containing formaldehyde as a fixative. Follow the instructions of the lysing solution manufacturer.
- Centrifuge tubes for 5 minutes at 300 g.
- Remove supernatant and resuspend pellet with 3 - 4 ml of PBS.
- Centrifuge tubes for 5 minutes at 300 g.
- Remove supernatant and resuspend pellet with 0.3 - 0.5 ml of PBS.
- Analyze samples immediately using flow cytometer or store samples at 2-8°C in the dark and analyze within 24 hours provided that cells were fixed.

Flow Cytometric Analysis

Analyze sample stained with CD56 PE using a flow cytometer. Visualize recorded data on the side-scatter (SSC) versus forward-scatter (FSC) plot. Set the gate for lymphocyte population as shown in figure 1. Then make a histogram of lymphocytes with PE intensity on the x-axis as shown in figure 2. Separate positive and negative populations using appropriate gates and calculate the percentage of CD56 positive lymphocytes. The region corresponding to the negative population should be set up using control cells which were stained by isotype control antibody.

Fig. 1: Delimitation of lymphocyte population





Fig. 2: Lymphocytes stained with CD56 PE reagent



10. Analytical performance

Specificity

The antibody MEM-188 reacts with a 180 kDa isoform of CD56 (NCAM) characteristic for leukocytes. CD56 is a transmembrane glycoprotein of immunoglobulin supergene family; the molecular weights of different isoforms ranging from 135 to 220 kDa. The monoclonal antibody MEM-188 was assigned to CD56 during the Human Leukocyte Differentiation Antigen workshop (HLDA6 WS Code: AS A055, NK NK26).

11. Clinical performance

Expected values

Results obtained in different laboratories may vary. Each laboratory should establish a normal range of cell subsets using its own test conditions. In our laboratory, the reagent CD56 PE was tested on 40 blood samples of healthy people. Obtained results are given in the table below.

Parameter	Mean (%)	SD	CV (%)
CD56+ lymphocytes	9.8	4.1	41.8

12. References

McCluggage WG et al. (2007) CD56 is a sensitive and diagnostically useful immunohistochemical marker of ovarian sex cord-stromal tumors. *Int J Gynecol Pathol.* 26: 322-327

Ohishi Y et al. (2007) CD56 expression in ovarian granulosa cell tumors, and its diagnostic utility and pitfalls. *Gynecol Oncol.* 107: 30-8

Zeromski J et al. (2005) Prevalence of CD56/NCAM molecule in nervous system, immune system and endocrine glands - accidental coincidence? *Endokrynol Pol.* 56: 78-82

Zeromski J et al. (2001) Significance of cell adhesion molecules, CD56/NCAM in particular, in human tumor growth and spreading. *Folia Histochem Cytobiol.* 39 Suppl. 2: 36-7

Leukocyte Typing VII., Mason D. et al. (Eds.), Oxford University Press (2002).

Leukocyte Typing VI., Kishimoto T. et al. (Eds.), Garland Publishing Inc. (1997).

13. Manufacturer

EXBIO Praha, a.s.
Nad Safinou II 341
25250 Vestec
Czech Republic

info@exbio.cz
technical@exbio.cz
orders@exbio.cz
www.exbio.cz

14. Trademarks

N/A

15. Revision History

- Version 1, ED7019_IFU_v1
Initial Release
- Version 2, ED7019_IFU_v2
Merging three language mutations into one document.
Intended use section was changed - "The reagent could be used in various antibody panels for multi-parameter flow cytometry analyses." - removed.
Staining protocol section was changed - step 5 - "(e.g. ADG-LYSE produced by AN DER GRUB, Cat. No. GAS-003)." - removed.
The address was changed: "Nad Safinou II 341."
- Version 3, ED7019_IFU_v3
Precautions section was changed - "Intended for professional use only." - removed. "Intended for In Vitro Diagnostic use in laboratories outside USA and Canada. This CE-IVD reagent is in conformity with the European In Vitro Diagnostic Medical Device Directive 98/79/EC." - added.
- Version 4, ED7019_IFU_v4
Reagent provided section was changed: text "stabilizing" added, "solution" - added and "0.2% (w/v) high-grade protease free Bovine Serum Albumin (BSA) as a stabilizing agent" - removed.
- Version 5, ED7019_IFU_v5
The company logo changed. IFU layout changed. "Keep away from sunlight." - added. Postal code changed: "25250 Vestec".

exbio

Monoclonal Antibody to CD56, PE conjugated (CD56 PE)

100 tests | Cat.No. ED7019



Instructions for Use

Version: ED7019_IFU_v5_EN

Date of Issue: 07-10-2020

EN

Symbols

The REF symbol consists of the letters 'REF' inside a rectangular border.	Catalogue number
The LOT symbol consists of the letters 'LOT' inside a rectangular border.	Batch code
The use-by date symbol is an hourglass icon.	Use-by date
The temperature limits symbol is a thermometer icon.	Temperature limits
The keep away from sunlight symbol is a sun with a slash through it.	Keep away from sunlight
The IVD symbol consists of the letters 'IVD' inside a rectangular border.	In vitro diagnostic medical device
The CE marking symbol consists of the letters 'C' and 'E' in a stylized font.	CE marking of conformity
The consult instructions for use symbol is an open book icon.	Consult instructions for use
The manufacturer symbol is a factory icon.	Manufacturer

The product is intended for In Vitro Diagnostic Use. In vivo diagnostic or therapeutic applications are strictly forbidden.

Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s..