CZ

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sHLA-G ELISA

Katalogové číslo:

RD194070100R

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16. REFERENCE

Reference k sHLA-G:

- Rebman V., Pfeiffer K., Passler M., Ferrone S., Maier S., Weiss E., Grosse-Wilde H.: Detection of soluble HLA-G molecules in plasma and amniotic fluid. *Tissue Antigens*. 1999; 53: 14-22
- Fuzzi B., Rizzo R., Criscuoli L., Noci I., Melchiorri L., Scarselli B., Bencini E., Menicucci A., Baricordi O.R.: HLA-G expression in early embryos is a fundamental prerequisite for the obtainment of pregnancy. *Eur J Immunol*. 2002; 32: 311-315
- Lila N., Amrein C., Guillemain R., Chevalier P., Latremouille C., Fabiani J.N., Dausset J., Carosella E.D., Carpentier A.: Human leukocyte antigen-G expression after heart transplantation is associated with a reduced incidence of rejection. *Circulation*. 2002; 105: 1949-1954
- Steinborn A., Rebmann V., Scharf A., Sohn C., Grosse-Wilde H.: Placental abruption is associated with decreased maternal plasma levels of soluble HLA-G. *J Clin Immunol*. 2003; 23: 307-314
- Le Bouteiller P., Pizzato N., Barakonyi A., Solier C.: HLA-G, pre-eclampsia, immunity and vascular events. *J Reprod Immunol*. 2003; 59: 219-234
- Sher G., Keskinetepe L., Nouriani M., Roussev R., Batzofin J.: Expression of sHLA-G in supernatants of individually cultured 46-h embryos: a potentially valuable indicator of 'embryo competency and IVF outcome. *Reprod Biomed Online*. 2004; 9: 74-78
- Noci I., Fuzzi B., Rizzo R., Melchiorri L., Criscuoli L., Dabizzi S., Biagiotti R., Pellegrini S., Menicucci A., Baricordi O.R.: Embryonic soluble HLA-G as a marker of developmental potential in embryos. *Hum Reprod*. 2005; 20: 138-146
- Yie S.M., Balakier H., Motamedi G., Librach C.L.: Secretion of human leukocyte antigen-G by human embryos is associated with a higher in vitro fertilization pregnancy rate. *Fertil Steril*. 2005; 83: 30-36
- Rebman V., Lemaoult J., Rouas-Friess N., Carosella E.D., Grosse-Wilde H.: Report of the Wet Workshop for Quantification of Soluble HLA-G in Essen, 2004. *Hum Immunol*. 2005; 66: 853-863
- Le Rond S., Azema C., Krawice-Radanne I., Durrbach A., Guettier C., Carosella E.D., Rouas-Freiss N.: Evidence to support the role of HLA-G5 in allograft acceptance through induction of immunosuppressive/ regulatory T cells. *J Immunol*. 2006; 176: 3266-3276
- Luque J., Torres M.I., Aumente M.D., Marin J., Garcia-Jurado G., Gonzalez R., Pascual D., Guerra N., Lopez-Rubio F., Alvarez-Lopez M.R., Arizon J.M., Pena J.: Soluble HLA-G in heart transplantation: their relationship to rejection episodes and immunosuppressive therapy. *Hum Immunol*. 2006; 67: 257-263
- Rudstein-Svetlicky N., Loewenthal R., Horejsi V., Gazit E.: HLA-G levels in serum and plasma. *Tissue Antigens* 2006; 67: 111-116
- Carosella E.D., Favier B., Rouas-Freiss N., Moreau P., Lemaoult J.: Beyond the increasing complexity of the immunomodulatory HLA-G molecule. *Blood*. 2008; 111(10):4862-70
- Apps R., Gardner L., Moffett A.: A critical look at HLA-G. *Trends Immunol*. 2008; 29(7):313-21
- Baricordi O.R., Stignani M., Melchiorri L., Rizzo R.: HLA-G and inflammatory diseases. *Inflamm Allergy Drug Targets*. 2008; 7(2):67-74

Reference k tomuto produktu:

- Wiendl H., Feger U., Mittelbronn M., Jack C., Schreiner B., Stadelmann C., Antel J., Brueck W., Meyermann R., Baror A., Kieseier B.C., Weller M.: Expression of the immune-tolerogenic major histocompatibility molecule HLA-G in multiple sclerosis: implications for CNS immunity. *Brain*. 2005; 128(Pt 11): 2689-2704
- Desai N., Filipovits J., Singh P., Beman M., Goldfarb J.: Assessment of soluble HLA-G levels in day 3 human embryo culture media using a new ELISA kit: Is there a relationship to embryo quality, development and implantation potential? *Fertility & Sterility*. 2005; 84, Suppl 1 (S281-S282)
- Steinborn A., Varkonyi T., Scharf A., Bahlmann F., Klee A., Sohn C.: Early detection of decreased soluble HLA-G levels in the maternal circulation predicts the occurrence of preeclampsia and intrauterine growth retardation during further course of pregnancy. *Am J Reprod Immunol* . 2007; 7(4):277-86
- Tahan F., Patirolu T. Plasma soluble human leukocyte antigen G levels in asthmatic children. *Int Arch Allergy Immunol* . 2006; 141(3):213-6 (2006)
- Desai N., Filipovits J., Goldfarb J.: Secretion of soluble HLA-G by day 3 human embryos associated with higher pregnancy and implantation rates: assay of culture media using a new ELISA kit. *Repro Biomed Online* . 2006; 13(2):272-7
- Coulam C.B., Roussev R.G.: Soluble HLA G concentrations in culture media from in vitro fertilized human embryos are predictive of implantation potential. *American Journal of Reproductive Immunology*. 2005; 53:279-309
- Wang F., Wen Z., Li H., Yang Z., Zhao X., Yao X.: Human leukocyte antigen-G is expressed by the eutopic and ectopic endometrium of adenomyosis. *Fertil Steril*. 2008; 90(5):1599-604
- Shaikly V.R., Morrison I.E., Taranissi M., Noble C.V., Withey A.D., Cherry R.J., Blois S.M. and Fernández N.: Analysis of HLA-G in maternal plasma, follicular fluid, and preimplantation embryos reveal an asymmetric pattern of expression. *J Immunol*. 2008; 180(6):4330-7
- Chen X.Y., Yan W.H., Lin A., Xu H.H., Zhang J.G., Wang X.X.: The 14 bp deletion polymorphisms in HLA-G gene play an important role in the expression of soluble HLA-G in plasma. *Tissue Antigens*. 2008; 72(4):335-41
- Ciprandi G., Colombo B.M., Contini P., Cagnati P., Pistorio A., Puppo F., Murdaca G.: Soluble HLA-G and HLA-A,-B,-C serum levels in patients with allergic rhinitis. *Allergy*. 2008; 63(10):1335-8
- Eidukaite A., Tamosiunas V.: Soluble HLA-G in the peritoneal fluid of women with endometriosis. *Fertil Steril*. 2008; 89(2):465-7
- Le Maux A., Noël G., Birebent B., Grosset J.M., Vu N., De Guibert S., Bernard M., Semana G., Amiot L: Soluble human leucocyte antigen-G molecules in peripheral blood haematopoietic stem cell transplantation: a specific role to prevent acute graft-versus-host disease and a link with regulatory T cells. *Clin Exp Immunol*. 2008; 152(1):50-6
- Giannopoulos K., Schmitt M., Kowal M., Własiuk P., Bojarska-Junak A., Roliński J., Dmoszyńska A.: The significance of soluble HLA-G plasma levels as well as messenger HLA-G for B-cell chronic lymphocytic leukemia (B-CLL). *Leuk Res*. 2008; 32(12):1815-9
- McCormick J., Whitley G.S., Le Bouteiller P., Cartwright J.E.: Soluble HLA-G regulates motility and invasion of the trophoblast-derived cell line SGHPL-4. *Hum Reprod*. 2009; 24(6):1339-45
- Kusanovic J.P., Romero R., Jodicke C., Mazaki-Tovi S., Vaisbuch E., Erez O., Mittal P., Gotsch F., Chaiworapongsa T., Edwin S.S., Pacora P., Hassan S.S.: Amniotic fluid soluble human leukocyte antigen-G in term and preterm parturition, and intra-amniotic infection/inflammation. *J Matern Fetal Neonatal Med*. 2009; 22(12):1151-66
- Roussev R.G., Acacio B., Parrett S., Romineb H., Coulam C.B.: Soluble HLA-G as potential candidate to help morphology grading score in embryo selection for IVF/ET. *Journal of Reproductive Immunology*. 2009; 81, 113–175
- Ivanova-Todorova E., Mourdjeva M., Kyurkchiev D., Bochev I., Stoyanova E., Dimitrov R., Timeva T., Yunakova M., Bukarev D., Shterev A., Tivchev P., Kyurkchiev S.: HLA-G

expression is up-regulated by progesterone in mesenchymal stem cells. Am J Reprod Immunol. 2009; 62(1):25-33

- Mach P., Blecharz P., Basta P., Marianowski P., Skret-Magierlo J., Kojs Z., Grabiec M., Wicherek L.: Differences in the soluble HLA-G blood serum concentration levels in patients with ovarian cancer and ovarian and deep endometriosis. Am J Reprod Immunol. 2010; 63(5):387-95
- Chen C.H., Liao H.T., Chen H.A., Liu C.H., Liang T.H., Wang C.T., Tsai C.Y., Chou C.T.: Human leukocyte antigen-G in ankylosing spondylitis and the response after tumour necrosis factor-alpha blocker therapy. Rheumatology (Oxford). 2010; 49(2):264-70

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