# EXCELLYSE® Easy 100 ml

## 10x concentrated lysing solution

REF ED7066

ENGLISH

### 1. Intended use

The EXCELLYSE® Easy lysing solution is intended for red blood cell lysis following antibody staining of human peripheral blood leukocytes.

The lysing solution is appropriate for use with EXBIO single colour monoclonal antibodies and KOMBITEST™ reagents and may be used in both lyse/wash and lyse/no wash protocol.

### 2. Introduction

Leukocyte analysis and detection in peripheral blood requires elimination of interfering cells, mainly erythrocyte. Ficoll density gradient method is usually used to separate leukocytes from whole blood. This method is rather time consuming and may lead to a loss of certain leukocyte subsets. Direct blood sample staining followed by red blood cell lysis therefore takes place in clinical laboratories as a fast and accurate method for whole blood flow cytometry analysis.

### 3. Reagent provided

The content of the vial (100 ml) is sufficient for 1000-2000 tests, respectively. The lysing solution is a 10× concentrated and **must be diluted 10× with deionized water prior to use (1 volume of concentrated solution with 9 volumes of deionized water).**

### 4. Storage

10× concentrated EXCELLYSE® Easy store at 2-25 °C. The prepared solution (1× concentrated) is stable for 1 month when stored at room temperature.

### 5. Precautions

- Intended for In Vitro Diagnostic use in laboratories outside USA and Canada. This CE-IVD reagent is in conformity with the European In Vitro Diagnostic Medical Device Directive 98/79/EC.
- Do not use after expiration date stamped on vial label.
- Do not use if any discoloration or precipitation occurs.
- Do not freeze.
- Avoid contamination of the reagent.
- Any non-performance of lysing protocol may produce false results.

- Blood samples are considered as potentially infectious and must be handled with care. Avoid all contact of the sample with the skin, eyes and mucosa.

**Warning:** The reagent contains formaldehyde, methanol and diethyleneglycol.

Danger

#### H-phrases

H302+H312+H332: Harmful if swallowed, in contact with skin or if inhaled.

H315: Causes skin irritation.

H317: May cause an allergic skin reaction.

H319: Causes serious eye irritation.

H335: May cause respiratory irritation.

H351: Suspected of causing cancer.

H371: May cause damage to organs.

H373: May cause damage to organs through prolonged or repeated exposure.

#### P-phrases

P270: Do not eat, drink or smoke when using this product.

P280: Wear protective gloves/protective clothing/eye protection/face protection.

P301 + P312: IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell.

P302 + P352: IF ON SKIN: Wash with plenty of soap and water.

P305 + P351 + P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P501: Dispose of contents/container to authorized facility for dangerous wastes.

### 6. Necessary material not supplied

Material necessary for collection of peripheral blood, test tubes for blood sample staining (e.g. 12 × 75 mm), automatic pipettes with disposable tips, vortex mixer, centrifuge, phosphate buffered saline (PBS), deionized water, appropriate fluorescent-dye-labeled primary/secondary antibody, flow cytometer.

## 7. Lyse/**no wash** lysing protocol

1. Collect peripheral blood in a sterile tube with an anticoagulant (e.g. Heparin, EDTA).
2. Follow instructions for whole blood antibody staining.
3. **Add 450 µl of diluted lysing solution per 50 µl of whole blood.** Mix the tube with a vortex mixer.
4. Incubate for about 5-10 minutes, until the blurry blood sample solution becomes clear.
5. Analyze the sample immediately using flow cytometer or store the sample at 2-8°C in the dark and analyze within 24 hours. No further cell fixation is required. See figures 1 and 2 for example data.

## 8. Lyse/**wash** lysing protocol

1. Collect peripheral blood in a sterile tube with an anticoagulant (e.g. Heparin, EDTA).
2. Follow manufacturer's instructions for whole blood antibody staining.
3. **Add 1000 µl of diluted lysing solution per 50 µl of whole blood.** Mix the tube with a vortex mixer.
4. Incubate for about 5-10 minutes, until the blurry blood sample solution becomes clear.
5. Centrifuge the tube for 5 minutes at 300 g.
6. Remove supernatant and resuspend the pellet with 0.2 – 0.5 ml of PBS.
7. Analyze the sample immediately using flow cytometer or store the sample at 2-8 °C in the dark and analyze within 24 hours. No further cell fixation is required. See figures 3 and 4 for example data.

## 9. Limitations

- Red blood cells from abnormal patients may be resistant to lysis using lysing solutions.
- Flow cytometer may produce false results if the device has not been aligned and maintained appropriately.
- Data may be incorrectly interpreted if fluorescent signals were compensated wrongly or if gates were positioned inaccurately.
- Blood samples from abnormal patients may exhibit abnormal values of positive cells.

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

**Fig. 1:** Peripheral blood leukocytes dot-plot from **lysed/non-washed** whole blood, analyzed on BD FACSCanto™ cytometer.

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
**Fig. 2:** Identification of CD45+ lymphocyte population in **lysed/non-washed** whole blood stained with anti-CD45 antibody, analyzed on BD FACSCanto™ cytometer.

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
**Fig. 3:** Peripheral blood leukocytes dot-plot from **lysed/washed** whole blood, analyzed on BD FACSCanto™ cytometer.

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**Fig. 4:** Identification of CD45+ lymphocyte population in **lysed/washed** whole blood stained with anti-CD45 antibody, analyzed on BD FACSCanto™ cytometer.

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