ENGLISH

1. Manufacturer

EXBIO Praha, a.s.
Nad Safinou II 341
252 42 Vestec, Czech Republic
Tel: +420 261 090 666
Fax: +420 261 090 660
E-mail: orders@exbio.cz
www.exbio.cz

2. Intended use

The SpermFlowEx[®] Kit is designed for analysis of following parameters in human semen using flow cytometry:

- **Sperm count**
- **Leukocyte count**
- **Sperm viability**
- **Sperm acrosome integrity**
- **Presence of intra-acrosomal protein in sperm**

The analysis should be completed by the examination of sperm motility and morphology using a light microscope in accordance to WHO recommendation [1].

3. Introduction

Lately, infertility of human population is a growing problem. Nearly 20 % of couples suffer from infertility, which is in 1/3 causes attributed to a male factor. Therefore, analysis of semen should play a role in basic non-invasive examination, which could verify or negate andrological cause of infertility.

Sperm examination is carried out mainly using the light microscope and obeys the WHO criteria from 2010 [1]. Such examination is subjective and depends on the personal experience of the examiner. When using the flow cytometry the analysis is balanced and objective, since the amount of measured sperms is much greater and detection of cells is provided via specific staining [2]. Main parameters having great impact on the ability of sperm to fertilize are: sperm count, sperm viability, acrosome integrity and presence of an intra-acrosomal protein (IAP) [3]. The acrosome contains digestive enzymes (including hyaluronidase and acrosin). These enzymes breakdown the outer membrane of the ovum called the zona pellucida, allowing the haploid nucleus of the sperm penetrate into the ovum. Presence of leukocytes in semen is a mark of an actual inflammation or a venereal disease.

4. Principle

Sperm count and leukocyte count

Measurement of the sperm count and leukocyte count is based on the addition of an internal standard (fluorescent beads with known concentration) to the semen sample. Detection of leukocytes in the sample is performed by staining with labeled antibody against human CD45 antigen (CD45 PE-Cy[™]5).

Sperm viability

Sperm viability is examined using propidium iodide which permeates through damaged membranes of dead cells and binds to their DNA.

Acrosome integrity

Measurement of acrosome integrity is based on the detection of an intra-acrosomal protein (IAP), which can be found inside the acrosome. If the sperm acrosome is intact, it is unable to detect IAP. Sperm with damaged membrane has IAP exposed and therefore accessible to the antibody against IAP, hence the protein is detected.

Presence of intra-acrosomal protein

After permeabilization of sperm membrane, intra-acrosomal protein is exposed to the antibody against intra-acrosomal protein (IAP) and thus is detected. In case, that sperm does not contain intra-acrosomal protein (IAP), the protein is not detected after permeabilization.

5. Precautions

- Intended for research use only.
- The flow cytometer should be calibrated on a routine basis using fluorescent microbeads to ensure stable sensitivity of detectors.
- Do not use reagents after expiration date stated on vial labels.
- Avoid prolonged exposure of the reagents to light.
- Avoid contamination of the reagents.
- Any non-performance of assay procedure may produce false results.
- Semen samples are considered as potentially infectious and must be handled with care.

Warning: The Permeabilizing Solution (ED7079-5) contains ethanol.

Warning

H-phrases

H226: Flammable liquid and vapour.

P-phrases

- P210: Keep away from heat/sparks/open flames/hot surfaces.
No smoking.
- P260: Do not breathe vapours/spray.
- P280: Wear protective gloves / eye protection / face protection.
- P301+P312: IF SWALLOWED: Call a POISON Center or doctor/physician if you feel unwell.
- P302+P352: IF ON SKIN: Wash with plenty of soap and water.
- P305+P351+P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
- P501 Dispose of contents/container to authorized facility for dangerous wastes.

6. Reagents provided

- ED7079-1 Intra-acrosomal protein FITC – 1 vial containing 0.5 ml of mouse monoclonal antibody against intra-acrosomal protein, FITC labeled.
- ED7079-2 CD45 PE-Cy™5 – 1 vial containing 0.25 ml of mouse monoclonal antibody against CD45, PE-Cy™5 labeled.
- ED7079-3 Fluorescent Count Standard – 1 vials containing 2 ml of fluorescent beads, 1×10^6 particles/ml
- ED7079-4 Propidium Iodide – 1 vials containing 0.25 ml of propidium iodide solution.
- ED7079-5 Permeabilizing Solution – 1 vials containing 25 ml of the reagent.

7. Storage

Store the SpermFlowEx® Kit at 2-8 °C. Expiration date is stated on vial labels and on the box.

8. Evidence of deterioration

In case of reagents deterioration or if data obtained show any performance alteration, please contact manufacturer using following e-mail address: technical@exbio.cz

9. Necessary material not supplied

- Suitable 5ml test tubes for blood staining (e.g. 12 × 75 mm)
- PBS buffer
- Automatic pipettes with disposable tips
- Vortex mixer
- Light microscope
- Centrifuge with rotor suitable for test tubes
- Flow cytometer - blue laser excitation at 488 nm and proper emission filters.

10. Specimen

Treat semen samples within 8 hours after ejaculation. Store semen samples at room temperature.

11. Assay procedure

Check the sperm count using a light microscope. Dilute the semen 1:1 or 1:9 with PBS according to the sperm count (the greater the count the greater the dilution).

Sperm and leukocyte count

1. Pipette 50 µl of diluted semen to a sample tube. Add 10 µl of **CD45 PE-Cy™5** reagent.
2. Mix the sample and incubate for 20 minutes at room temperature.
3. Add 50 µl of **Fluorescent Count Standard** (first drop well mixed standard in an empty microtube and then pipette the exact volume of the standard to the sperm sample).
4. Add 0.5 ml of PBS, mix the sample and analyze using a flow cytometer.

Sperm viability

1. Pipette 50 µl of diluted semen to a sample tube.
2. Add 10 µl of **Propidium Iodide** solution.
3. Mix the sample and incubate for 20 minutes at room temperature.
4. Add 0.5 ml of PBS, mix the sample and analyze using a flow cytometer.

Acrosome integrity

1. Pipette 50 µl of diluted semen to a sample tube.
2. Add 1 ml of PBS.
3. Mix well and centrifuge for 5 minutes at 150 g.
4. Wash the pellet twice (repeat the 2nd and 3rd step).
5. Add 10 µl of **Intra-acrosomal protein FITC** antibody to the pellet.
6. Mix the sample and incubate for 20 minutes at room temperature.
7. Add 0.5 ml of PBS, mix well and analyze using a flow cytometer.

Presence of intra-acrosomal protein

1. Pipette 50 µl of diluted semen to a sample tube.
2. Add 1 ml of **Permeabilizing Solution** which **was allowed to warm up to room temperature** (20-25 °C).
3. Mix the sample and incubate for 15 minutes at room temperature.
4. Centrifuge for 5 minutes at 150 g.
5. Add 1 ml of PBS to the pellet, mix and centrifuge for 5 minutes at 150 g.
6. Wash the pellet twice (repeat the 5th step).
7. Add 10 µl of **Intra-acrosomal protein FITC** antibody to the pellet.
8. Mix well and incubate for 20 minutes at room temperature.
9. Add 0.5 ml of PBS, mix well and analyze using a flow cytometer.

12. Flow cytometry

Analyze stained samples using flow cytometer equipped with a 488 nm excitation laser and appropriate filter set-up.

Fluorescent Count Standard is detected in FITC fluorescence detector as well as monoclonal antibody against intra-acrosomal protein. Fluorescence of propidium iodide and PE-Cy™5 label is detected in PC5 or PerCP fluorescence detector.

13. Data analysis

Sperm and leukocyte count

Visualize measured data as a dot-plot, where FITC fluorescence intensity is on the X-axis and forward-scatter (FSC) is on the Y-axis (Fig.1). Set gates as shown in figure 1: gate **A** for all cells, gate **B** for sperms and gate **C** for Fluorescent Count Standard (fluorescent beads). Then visualize events from the gate **A** as a dot-plot (Fig. 2), where the X-axis represents fluorescence intensity of PE-Cy™5 dye and the Y-axis represents side-scatter (SSC). On the contrary to the sperms, leukocytes are CD45 positive. Separate the bright population of leukocytes using gate **D** as shown in figure 2. Calculate sperm and leukocyte count using following formulas:

$$S = \frac{\text{gateB}}{\text{gateC}} \times \text{dilution} \quad [10^6 \text{ sperm/ml}]$$

$$L = \frac{\text{gateD}}{\text{gateC}} \times \text{dilution} \quad [10^6 \text{ leukocytes/ml}]$$

S represents sperm count in ejaculate
L represents leukocyte count in ejaculate
Gate B represents event count in gate B
Gate C represents event count in gate C
Gate D represents event count in gate D

Sperm viability

Visualize measured data as a forward-scatter (FSC) versus side-scatter (SSC) dot-plot. Set the gate around sperm cells as shown in figure 3. Then visualize sperm cells in a histogram (Fig. 4), where the X-axis represents fluorescence intensity of propidium iodide in PE-Cy5 or PerCP channel. Separate positive and negative sperm cells using appropriate gate. Negative population represents viable sperms, while PI-positive population represents non-viable sperms. Viability of sperms is represented by the percentage of viable sperms from all sperm cells.

Acrosome integrity and presence of intra-acrosomal protein

Visualize permeabilized and non-permeabilized samples stained using antibody against intra-acrosomal protein as a forward-scatter (FSC) versus side-scatter (SSC) dot-plot and set the gate around sperm events as shown in figure 3. Then visualize sperm events in a histogram (see Fig. 5-7), where the X-axis represents fluorescence intensity in FITC fluorescence detector. Separate positive and negative sperms using appropriate gates.

Figure 5 shows non-permeabilized non-pathological sperms events with low content (< 30 %) of positive sperms. Negative ones represent sperms with intact acrosome and positive ones represent sperms with damaged acrosome.

Figure 6 shows permeabilized non-pathological sperms with dominant content of positive sperm events. Positive sperm events indicate the presence of the intra-acrosomal protein.

Figure 7 shows non-permeabilized pathological sperms with high content (> 30 %) of positive sperm events with damaged acrosome.

14. Expected values

Normal values according to WHO [1]

Sperm count in semen > 15×10⁶/ml
Leukocyte count in semen < 1×10⁶/ml
Sperm viability > 58 %

Normal acrosomal values

Acrosome integrity of sperms < 30 % pathological
Presence of intra-acrosomal protein > 90 % sperms

15. Limitations

- Composition of components was optimized to offer the best specific signal/non-specific signal ratio. Therefore, it is important to adhere to the reagent volume/sample volume ratio in every test.
- Flow cytometer may produce false results if the device has not been aligned and maintained appropriately.
- In case of low sperm counts (less than 10 x 10⁶ /ml) we recommend to check sperm count using microscope. In such samples the results might be affected due to relatively higher content of cell debris in the gated region.
- The monoclonal antibody against Intra-acrosomal protein cannot be used for specific identification of human sperms due to the antibody binding to cell debris or potential cross reactivity with other isoforms of target antigen (sperm-specific glyceraldehyde phosphate dehydrogenase) present in other cell types.

16. Trademarks

Cy™ and CyDye™ are registered trademarks of GE Healthcare.

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18. References / Literatura / Literatúra


[1] WHO laboratory manual for the Examination and processing of human semen. World Health Organization, 5th edition, 2010

[2] Gilan L, Evans G, Maxwell WM (2005) Flow cytometric evaluation of sperm parameters in relation to fertility potential. Theriogenology. 63: 445-57

[3] Peknicova J, Chladek D, Hozak P (2005) Monoclonal antibodies against sperm intra-acrosomal antigens as markers for male infertility diagnostics and estimation of spermatogenesis. Am J Reprod Immunol. 53(1): 42-9

[4] Leukocyte Typing IV., Knapp W. et al. (Eds.), Oxford University Press (1989).

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