

Monoclonal Antibody to CD45, PerCP conjugated (CD45 PerCP)

Cat.No. ED7029

1. Intended purpose

The reagent CD45 PerCP permits identification and enumeration of cell populations expressing human CD45 antigen in whole blood using flow cytometry.

2. Test principle

This test is based on specific binding of monoclonal antibody to the antigenic determinant expressed on the surface of leukocytes. The monoclonal antibody is labeled with fluorochrome which is excited via laser beam from a flow cytometer during analysis. Subsequent emission of light from fluorochromes of each cell is collected and analyzed by a flow cytometer. The fluorescence intensity differences enable separation of cell subsets based on expression of analyzed antigen. Specific staining of blood cells is performed by incubation of blood samples with the reagent followed by a lysis of red blood cells. Afterwards, unaffected leukocytes are subjected to analysis by a flow cytometer.

3. Reagents provided

The reagent contains mouse monoclonal antibody against human CD45 antigen (clone MEM-28) which was purified by affinity chromatography and labeled with Peridinin-chlorophyll-protein complex (PerCP). The labeled antibody is diluted at optimum concentration in stabilizing phosphate buffered saline (PBS) solution containing 15mM sodium azide. The content of a vial (2 ml) is sufficient for 100 tests.

Product specification

Content	100 tests, 2 ml
Usage	20 µl per test
Specificity	Human CD45
Clone	MEM-28
Isotype	Mouse IgG1
Fluorochrome	PerCP
λ excitation	488 nm
Emission maximum	670 nm

4. Materials required but not provided

Test tubes for blood staining (e.g. 12 × 75 mm)
Commercial lysing solution
Phosphate buffered saline (PBS)
Isotype control antibody (mouse IgG1 PerCP)

5. Equipment required

Automatic pipettes with disposable tips
Vortex mixer
Centrifuge
Flow cytometer with excitation laser 488 nm and proper filters

6. Storage and handling

Store the vial at 2-8 °C. Keep away from sunlight. Do not freeze. Do not aliquote. Expiration date is stated on a vial label and on outer packaging.

7. Warnings, precautions and limitations of use

- Intended for In Vitro Diagnostic use in laboratories outside USA and Canada. This CE-IVD kit is in conformity with the European Directive 98/79/EC.
- Do not use reagent after expiration date.
- Avoid reagents contamination.
- Avoid prolonged exposure to light.
- The content of the vial must not freeze.
- Any non-performance of staining protocol may produce false results.
- The reagent contains sodium azide (NaN₃) which is highly toxic in pure form. However, the concentration in the reagent (15mM) is not considered as hazardous. When disposing the reagent, flush the sink with a large volume of water.
- Blood samples are considered as potentially infectious and must be handled with care. Avoid all contact of the sample with the skin, eyes and mucosa.
- In case of hyperleukocytose sample, it is recommended to dilute blood sample with PBS to obtain leukocyte density approximately 5 × 10⁶ leukocytes/ml.
- Blood samples from abnormal patients may exhibit abnormal values of positive cells.
- Data may be incorrectly interpreted if fluorescent signals were compensated

wrongly or if gates were positioned inaccurately.

- Flow cytometer may produce false results if the device has not been aligned and maintained appropriately.
- Red blood cells from abnormal patients may be resistant to lysis using lysing solutions.
- Blood samples should be stained and analyzed within 24 hours from the blood collection.

8. Specimen

Use the peripheral human blood in a sterile tube with an anticoagulant (Heparin or EDTA). Blood must be stored at room temperature.

9. Procedure


Staining protocol

- Add 20 µl of CD45 PerCP reagent to a test tube, and the necessary amount of isotype control to a control tube.
- Add 100 µl of blood sample to each tube. Vortex the tubes.
- Incubate tubes for 20-30 minutes at room temperature in the dark.
- Perform lysis of red cells using lysing solution. It is recommended to use a commercial lysing solution containing formaldehyde as a fixative. Follow the instructions of the lysing solution manufacturer.
- Centrifuge tubes for 5 minutes at 300 g.
- Remove supernatant and resuspend pellet with 3-4 ml of PBS.
- Centrifuge tubes for 5 minutes at 300 g.
- Remove supernatant and resuspend pellet with 0.3 - 0.5 ml of PBS.
- Analyze samples immediately using flow cytometer or store samples at 2-8 °C in the dark and analyze within 24 hours provided that cells were fixed.

Flow Cytometric Analysis

Analyze the sample stained with CD45 PerCP using a flow cytometer. Visualize recorded data using appropriate plot such as side-scatter (SSC) versus PerCP intensity as shown in figure 1. All leukocytes are bright (CD45+), nonleukocytes (debris, erythrocytes, platelets, etc.) are dim (CD45-). Set suitable gates for analysis.

Fig. 1: Leukocytes stained with CD45 PerCP reagent



10. Analytical performance

Specificity

The antibody MEM-28 reacts with all alternative forms of human CD45 phosphotyrosine phosphatase (Leukocyte Common Antigen), a 180-220 kDa single chain type I transmembrane protein expressed at high level on all cells of hematopoietic origin, except erythrocytes and platelets. The monoclonal antibody MEM-28 was assigned to CD45 during the Human Leukocyte Differentiation Antigen workshop (HLDA3 WS Code: NL 833a).

11. Clinical performance

N/A

12. References

- Guttinger M et al. (1992) CD45 phosphotyrosine phosphatase and p56lck protein tyrosine kinase: a functional complex crucial in T cell signal transduction. *Int Immunol.* 4: 1325-30.
- Stover DR et al. (1991) Protein-tyrosine-phosphatase CD45 is phosphorylated transiently on tyrosine upon activation of Jurkat T cells. *Proc Natl Acad Sci USA* 88: 7704-7707.
- Taetle R et al. (1991) Regulation of CD45 expression in human leukemia cells. *Leukemia* 5: 309-314.
- Nakano A et al. (1990) Expression of leukocyte common antigen (cd45) on various human leukemia/lymphoma cell lines. *Acta Pathol Jpn.* 40: 107-15.
- Yamada A et al. (1990) Effect of activation of protein kinase C on CD45 isoform expression and CD45 protein tyrosine phosphatase activity in T cells. *Eur J Immunol.* 20: 1655-60.
- Bazil V et al. (1989) Sialic acid-dependent epitopes of CD45 molecules of restricted cellular

expression. *Immunogenetics* 29: 202-5.

Horejsi V et al. (1988) Monoclonal antibodies against human leucocyte antigens. II. Antibodies against CD45 (T200), CD3 (T3), CD43, CD10 (CALLA), transferrin receptor (T9), a novel broadly expressed 18-kDa antigen (MEM-43) and a novel antigen of restricted expression (MEM-74). *Folia Biol (Praha).* 34: 23-34.

Leukocyte Typing IV., Knapp W. et al. (Eds.), Oxford University Press (1989).

Leukocyte Typing III., McMichael A. J. et al (Eds.), Oxford University Press (1987).

13. Manufacturer

EXBIO Praha, a.s.
Nad Safinou II 341
25250 Vestec
Czech Republic

info@exbio.cz
technical@exbio.cz
orders@exbio.cz
www.exbio.cz

14. Trademarks

N/A

15. Revision History

- Version 1, ED7029_IFU_v1
Initial Release
- Version 2, ED7029_IFU_v2
The address was changed: "Nad Safinou II 341"
- Version 3, ED7029_IFU_v3
Precautions section was changed - "Intended for professional use only." - removed. "Intended for In Vitro Diagnostic use in laboratories outside USA and Canada. This CE-IVD reagent is in conformity with the European In Vitro Diagnostic Medical Device Directive 98/79/EC." - added.
- Version 4, ED7029_IFU_v4
Precautions section was changed - stabilizing added, solution - added and "0.2% (w/v) high-grade protease free Bovine Serum Albumin (BSA) as a stabilizing agent" - removed.
- Version 5, ED7029_IFU_v5
The company logo changed. IFU layout changed. "Keep away from sunlight." - added. Postal code changed: "25250 Vestec"

exbio

Monoclonal Antibody to CD45, PerCP conjugated (CD45 PerCP)

100 tests | Cat.No. ED7029



Instructions for Use

Version: ED7029_IFU_v5_EN

Date of Issue: 27-02-2020

EN

Symbols

The REF symbol consists of the letters 'REF' in a bold, sans-serif font, enclosed within a rectangular border.	Catalogue number
The LOT symbol consists of the letters 'LOT' in a bold, sans-serif font, enclosed within a rectangular border.	Batch code
The use-by date symbol is an icon of an hourglass.	Use-by date
The temperature limits symbol is an icon of a thermometer.	Temperature limits
The keep away from sunlight symbol is an icon of a sun with rays.	Keep away from sunlight
The IVD symbol consists of the letters 'IVD' in a bold, sans-serif font, enclosed within a rectangular border.	In vitro diagnostic medical device
The CE marking symbol consists of the letters 'C' and 'E' in a bold, sans-serif font, enclosed within a partial rectangular border.	CE marking of conformity
The consult instructions for use symbol is an icon of an open book.	Consult instructions for use
The manufacturer symbol is an icon of a factory building.	Manufacturer

The product is intended for In Vitro Diagnostic Use. In vivo diagnostic or therapeutic applications are strictly forbidden.

Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s..