KOMBITEST™ CD3 FITC / CD4 PE-Cy™7 / CD8 APC-Cy™7 / CD16+CD56 PE / CD19 APC / CD45 PerCP-Cy™5.5

Cat. No. / Kat. č.: ED7075

ENGLISH

1. Manufacturer

EXBIO Praha, a.s.
Nad Safinou II 341
252 42 Vestec, Czech Republic
Tel: +420 261 090 666
Fax: +420 261 090 660
E-mail: orders@exbio.cz
www.exbio.cz

2. Specification

Content	50 tests, 1 ml						
Usage	20 µl per test						
Specificity	CD3	CD4	CD8	CD16	CD56	CD19	CD45
Clone	UCHT1	MEM-241	MEM-31	3G8	LT56	LT19	MEM-28
Isotype (mouse)	IgG1	IgG1	IgG2s	IgG1	IgG2a	IgG1	IgG1
Fluorochrome	FITC	PE-Cy™7	APC-Cy™7	PE		APC	PerCP-Cy™5.5
λ. excitation	488 nm	488 nm	633 nm	488 nm		633 nm	488 nm
Emission maximum	525 nm	780 nm	780 nm	575 nm		660 nm	692 nm

3. Intended use

The KOMBITEST™ CD3 FITC / CD4 PE-Cy™7 / CD8 APC-Cy™7 / CD16+CD56 PE / CD19 APC / CD45 PerCP-Cy™5.5 is designed for identification and enumeration of mature human T, B, and NK lymphocytes as well as helper/inducer (CD3+CD4+) and suppressor/cytotoxic (CD3+CD8+) T lymphocyte subsets in erythrocyte-lyzed whole blood using Flow Cytometry.

4. Principle

This test is based on the specific binding of monoclonal antibodies to the antigenic determinants expressed on the surface of leukocytes. The monoclonal antibodies are labeled with different fluorochromes which are excited via laser beam from a flow cytometer during analysis. Subsequent emissions of light from the fluorochromes of each cell are collected and analyzed by a flow cytometer. The fluorescence intensity differences enable the separation of cell subsets based on the expression of analyzed antigens.

The specific staining of blood cells is performed by the incubation of blood samples with the reagent followed by a lysis of red blood cells. Afterwards, unaffected leukocytes are subjected to analysis by a flow cytometer.

5. Reagent provided

The reagent contains premixed combination of mouse monoclonal antibody against human CD3 antigen (clone UCHT1) labeled with Fluorescein isothiocyanate (FITC), mouse monoclonal antibody against human CD4 antigen (clone MEM-241) labeled with tandem dye PE-Cy™7, mouse monoclonal antibody against human CD8 antigen (clone MEM-31) labeled with tandem dye APC-Cy™7, mouse monoclonal antibody against human CD16 antigen (clone 3G8) labeled with R-phycoerythrin (PE), mouse monoclonal antibody against human CD56

antigen (clone LT56) labeled with R-phycoerythrin (PE), mouse monoclonal antibody against human CD19 antigen (clone LT19) labeled with Allophycocyanin (APC), and mouse monoclonal antibody against human CD45 antigen (clone MEM-28) labeled with tandem dye PerCP-Cy™5.5. Labeled antibodies are diluted at optimum concentration in stabilizing phosphate buffered saline (PBS) solution containing 15mM sodium azide. The content of a vial (1 ml) is sufficient for 50 tests.

6. Storage

Store vial in the dark at 2-8°C. Do not freeze.

7. Evidence of deterioration

In case of reagent deterioration or if obtained data show any performance alteration, please contact manufacturer using following e-mail address: technical@exbio.cz

8. Precautions

- Intended for In Vitro Diagnostic use in laboratories outside USA and Canada. This CE-IVD reagent is in conformity with the European In Vitro Diagnostic Medical Device Directive 98/79/EC.
- Do not use after expiration date stamped on vial label.
- Avoid prolonged exposure to light.
- The content of the vial must not freeze.
- Avoid contamination of the reagent.
- The flow cytometer should be calibrated on a routine basis using fluorescent microbeads to ensure stable sensitivity of detectors.
- Any non-performance of the staining protocol may produce false results.
- The reagent contains sodium azide (NaN₃) which is highly toxic in pure form. However, the concentration in the reagent (15mM) is not considered as hazardous. When disposing the reagent, flush the sink with a large volume of water.
- Blood samples are considered as potentially infectious and must be handled with care. Avoid all contact of the sample with the skin, eyes and mucosa.

9. Necessary material not supplied

Test tubes for staining of blood samples (e.g. 12 × 75 mm), automatic pipettes with disposable tips, vortex mixer, EXCELLYSE Easy lysing solution (Cat.No. ED7066), flow cytometer with two lasers (488 nm and 633 or 635 nm).

10. Blood specimen

The peripheral blood in a sterile tube with an anticoagulant (Heparin or EDTA) must be stored at room temperature and stained within 48 hours of drawing.

11. Staining protocol

- Add 20 µl of KOMBITEST™ CD3 FITC / CD4 PE-Cy™7 / CD8 APC-Cy™7 / CD16+CD56 PE / CD19 APC / CD45 PerCP-Cy™5.5 reagent to a test tube.
- Add 50 µl of blood sample to the tube. Vortex the tube.
- Incubate the tube for 15-20 minutes at room temperature in the dark.
- Perform lysis of red cells using EXCELLYSE Easy lysing solution (Cat.No. ED7066) or any other commercial lysing

solution containing formaldehyde as a fixative. Follow the instructions of the lysing solution manufacturer.

- Analyze the sample immediately using a flow cytometer or store sample at 2-8°C in the dark and analyze within 24 hours provided that cells were fixed.

12. Flow cytometry

Analyze stained samples using a flow cytometer equipped with two excitation lasers (488 nm and 633 or 635 nm) and proper filters. Compensate fluorescent signals prior to or after data acquisition.

13. Data analysis

Visualize compensated data on the side-scatter (SSC) versus CD45 PerCP-Cy[™]5.5 plot. Set the gate for CD45+ lymphocyte population as shown in figure 1 (chapter 18).

Visualize CD45+ lymphocytes on the side-scatter (SSC) versus CD3 FITC plot as shown in figure 2 (chapter 18). Separate CD3+ and CD3- lymphocytes using appropriate gates. Calculate the percentage of T lymphocytes (CD3+).

Then plot the gated CD3- lymphocytes as CD19 APC versus CD16+CD56 PE as shown in figure 3 (chapter 18). Set appropriate gates and calculate the percentage of B lymphocytes (CD16-CD56-CD19+) and natural killer (NK) lymphocytes (CD16+CD56+CD19-). Please note, calculate the frequency of grandparent to get relative numbers in lymphocyte population.

Finally, plot T lymphocytes (CD3+) as CD4 PE-Cy[™]7 versus CD8 APC-Cy[™]7 as shown in figure 4 (chapter 18). Set appropriate gates and calculate the percentage of helper/inducer T lymphocytes (CD4+CD8-) and suppressor/cytotoxic T lymphocytes (CD4-CD8+). Please note, calculate the frequency of grandparent to get relative numbers in lymphocyte population.

14. Expected values

Results obtained in different laboratories may vary. Each laboratory should establish a normal range of cell subsets using its own test conditions. Results obtained in our laboratory are given in the table below.

Lymphocyte Subset	Unit	n	Mean	95% Range
CD3+	%	108	71	52-83
	cells/ μ l	54	1480	886-2341
CD3+CD8+	%	108	25	13-46
	cells/ μ l	54	503	208-1031
CD3+CD4+	%	108	44	25-64
	cells/ μ l	54	939	473-1524
CD16+CD56+	%	108	15	6-34
	cells/ μ l	54	319	121-536
CD19+	%	108	13	5-24
	cells/ μ l	54	242	87-445

15. Performance characteristics

Specificity

The antibody UCHT1 recognizes the CD3 antigen of the TCR/CD3 complex on mature human T cells. The UCHT1 antibody reacts with the epsilon chain of the CD3 complex.

HLDA I; WS Code T 3

HLDA III; WS Code T 126

HLDA III; WS Code T 471

HLDA VI; WS Code T 6T-CD3.1

The antibody MEM-241 recognizes CD4 co-receptor, a 55 kDa transmembrane glycoprotein of immunoglobulin family expressed on subsets of T lymphocytes (such as "helper" T cells, CD4+ regulatory T cells or CD4+CD8+ double-positive T cells) and also on monocytes, tissue macrophages and granulocytes.

HLDA VIII (HCDM); WS Code M241

The antibody MEM-31 recognizes a conformationally-dependent epitope of CD8, a cell surface glycoprotein that mediates efficient cell-cell interactions within the immune system. CD8 is a disulfide-linked dimer (each monomer approx. 32-34 kDa) and exists as a CD8 α /

homodimer on subsets of memory T cells, intraepithelial lymphocytes, NK cells and dendritic cells, or as a CD8 α / β heterodimer on majority of MHC I-restricted cytotoxic T cells and thymocytes.

HLDA III; WS Code T 575

The antibody 3G8 recognizes CD16, a low affinity receptor for aggregated IgG (Fc γ RIII antigen). CD16 exists in two different isoforms: CD16a (Fc γ RIIIA; 50-65 kDa; expressed on NK-cells, monocytes and macrophages) and CD16b (Fc γ RIIIB; 48 kDa; mainly expressed on neutrophils).

HLDA V; WS Code NK80

The antibody LT56 reacts with a 180 kDa isoform of CD56 (NCAM) characteristic for leukocytes. CD56 is a single transmembrane glycoprotein of immunoglobulin supergene family; the molecular weights of different isoforms ranging from 135 to 220 kDa.

The antibody LT19 recognizes CD19 (B4), a 95 kDa type I transmembrane glycoprotein of immunoglobulin superfamily, expressed on B lymphocytes and follicular dendritic cells; it is lost on plasma cells.

The antibody MEM-28 reacts with all alternative forms of human CD45 phosphotyrosine phosphatase (Leukocyte Common Antigen), a 180-220 kDa single chain type I transmembrane protein expressed at high level on all cells of hematopoietic origin, except from erythrocytes and platelets.

HLDA III; WS Code NL 833a

Accuracy

The accuracy of the method was studied by the comparison of KOMBITEST[™] CD3 FITC / CD4 PE-Cy[™]7 / CD8 APC-Cy[™]7 / CD16+CD56 PE / CD19 APC / CD45 PerCP-Cy[™]5.5 with competitor's product in parallel staining of blood samples. Samples were analyzed using BD FACSCanto[™] Flow Cytometer. The regression parameters are given in the table below.

Lymphocyte Subset	Unit	n	Slope	Intercept	r ²	Range
CD3+	%	109	0.987	2.26	0.949	48-86
	cells/ μ l	54	0.967	109	0.920	768-3081
CD3+CD8+	%	109	1.029	0.62	0.960	7.9-42
	cells/ μ l	54	1.064	9.32	0.950	203-1296
CD3+CD4+	%	109	0.956	2.87	0.954	25-66
	cells/ μ l	54	0.904	114	0.899	439-1813
CD16+CD56+	%	109	0.981	0.06	0.977	2.8-45
	cells/ μ l	54	1.047	-8.49	0.963	104-1099
CD19+	%	109	0.984	-0.26	0.951	2.1-27
	cells/ μ l	54	1.001	2.15	0.925	72-441

Linearity

The linearity of the method was verified on 10 serial dilutions of leukocyte-enriched blood sample (buffy coat). Cell samples were stained by KOMBITEST[™] CD3 FITC / CD4 PE-Cy[™]7 / CD8 APC-Cy[™]7 / CD16+CD56 PE / CD19 APC / CD45 PerCP-Cy[™]5.5 in triplicates. Measured data were observed to be linear across the tested range 50-25000 lymphocytes/ μ l. Cell subsets were in following ranges.

Lymphocyte Subset	Range (cells/ μ l)
CD3+	37-19028
CD3+CD8+	12-6156
CD3+CD4+	23-11996
CD16+CD56+	9-4804
CD19+	3-1646

Repeatability

The repeatability of the assay was measured on one blood sample in twelve tubes in parallel. Coefficients of variation (CV) are given in the table below.

Lymphocyte Subset	Unit	n	Average	SD	CV
CD3+	%	12	71.1	0.80	1.13
	cells/ μ l	12	2007	209	10.4
CD3+CD8+	%	12	27.0	0.97	3.59
	cells/ μ l	12	762	94	12.2
CD3+CD4+	%	12	42.6	0.91	2.14
	cells/ μ l	12	1203	122	10.1
CD16+CD56+	%	12	16.3	0.80	4.92
	cells/ μ l	12	461	62	13.4
CD19+	%	12	11.0	0.48	4.38
	cells/ μ l	12	310	27	8.67

Reproducibility

The reproducibility of the assay was measured on stabilized blood sample (Immuno-Troll™ Cells, Beckman-Coulter) under the same conditions for one month. Coefficients of variation (CV) are given in the table below.

Lymphocyte Subset	Unit	n	Average	SD	CV
CD3+	%	11	70.0	3.98	5.69
	cells/ μ l	11	865	63.4	7.33
CD3+CD8+	%	11	24.4	0.76	3.14
	cells/ μ l	11	301	14.7	4.87
CD3+CD4+	%	11	50.2	0.49	0.98
	cells/ μ l	11	620	20.8	3.35
CD16+CD56+	%	11	11.0	0.62	5.61
	cells/ μ l	11	135	6.63	4.89
CD19+	%	11	12.3	0.69	5.57
	cells/ μ l	11	152	9.58	6.31

16. Limitations

- Concentrations of labeled antibodies in this reagent were optimized to offer the best specific signal/non-specific signal ratio. Therefore, it is important to adhere to the reagent volume/sample volume ratio in every test.
- Flow cytometer may produce false results if the device has not been aligned and maintained appropriately.
- Data may be incorrectly interpreted if fluorescent signals were compensated wrongly or if gates were positioned inaccurately.
- Blood samples from abnormal patients may exhibit abnormal values of positive cells.
- In case of a hyperleukocytose sample, it is recommended to dilute the blood sample with PBS to obtain leukocyte density approximately 5×10^6 leukocytes/ml.

17. Trademarks

Cy™ and CyDye™ are registered trademarks of GE Healthcare. KOMBITEST™ is registered trademark of EXBIO Praha, a.s. BD FACSCanto™ is registered trademark of Becton, Dickinson and Company.

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