EXCELLYSE Live 100 ml

10× concentrated lysing solution

REF ED7068

ENGLISH

1. Intended use

The EXCELLYSE Live lysing solution is intended for red blood cell lysis following antibody staining of leukocytes in human peripheral whole blood. This solution is appropriate for use when viable leukocytes are required after red blood cell lysis.

The lysing solution is optimized for use with EXBIO single colour monoclonal antibodies and KOMBITEST™ reagents and may be used in both lyse/wash and lyse/no wash protocol.

2. Introduction

Leukocyte analysis and detection in peripheral blood requires elimination of interfering cells, mainly erythrocytes. So far, the Ficoll density gradient method was used to separate leukocytes from whole blood. This method is rather time consuming and may lead to a loss of certain cell subsets. Direct blood sample staining followed by red blood cell lysis therefore takes place in clinical laboratories as a fast and easy method for whole blood flow cytometry analysis.

The EXCELLYSE Live lysing solution allows for red blood cell lysis and is free of fixation component keeping leukocyte membrane intact and the cells viable.

3. Reagent provided

The content of the vial (100 ml) is sufficient for 500 tests, respectively. The lysing solution is a 10× concentrated and **must be diluted 10× with deionized water prior to use (1 volume of concentrated solution with 9 volumes of deionized water)**. Always prepare a fresh 1× solution before lysis.

4. Storage

Store the 10× concentrated EXCELLYSE Live at 2-8 °C. Freshly prepared solution (1× concentrated) is stable for 1 day when stored at room temperature.

5. Precautions

- For research use only.
- Do not use after expiration date stated on vial label.
- Do not use if any discoloration or precipitation occurs.
- Do not freeze.
- Avoid contamination of the reagent.
- Any non-performance of lysing protocol may produce false results.

- Blood samples are considered as potentially infectious and must be handled with caution. Avoid any contact of the sample with the skin, eyes and mucosa.

Warning: The reagent contains ammonium chloride which may produce a characteristic non-pleasant odor of ammonia when opened after a long-term storage.

6. Necessary material not supplied

Material necessary for collection of peripheral blood, test tubes for blood sample staining (e.g. 12 × 75 mm), automatic pipettes with disposable tips, vortex mixer, centrifuge, phosphate buffered saline (PBS), deionized water, appropriate fluorescent-dye-labeled primary/secondary antibody, flow cytometer.

7. Lyse/no wash blood lysing protocol

1. Collect peripheral blood in a sterile tube with an anticoagulant (e.g. Heparin, EDTA).
2. Follow instructions for the whole blood antibody staining.
3. **Add 2.0 ml of the diluted lysing solution per 100 µl of whole blood.** Mix the tube with a vortex mixer.
4. Incubate for 5-10 minutes, until the blurry blood sample solution becomes clear.
5. Analyze the sample immediately using flow cytometer or store the sample at 2-8 °C in the dark and analyze within 24 hours. See figures 1 and 2 for example data.

8. Lyse/wash blood lysing protocol

1. Collect peripheral blood in a sterile tube with an anticoagulant (e.g. Heparin, EDTA).
2. Follow manufacturer's instructions for the whole blood antibody staining.
3. **Add 2.0 ml of the diluted lysing solution per 100 µl of whole blood.** Mix the tube with a vortex mixer.
4. Incubate for 5-10 minutes, until the blurry blood sample solution becomes clear.
5. Centrifuge the tube for 5 minutes at 300×g.
6. Remove supernatant and resuspend the pellet with 0.2 – 0.5 ml of PBS.
7. Analyze the sample immediately using flow cytometer or store the sample at 2-8 °C in the dark and analyze within 2 hours.
8. Fix the cells if samples are required to be stored for more than 2 hours. Non-fixed cells may change their size. Note that cells lose their viability after fixation. See figures 3 and 4 for example data.

9. Limitations

- Red blood cells from some patients may be resistant to lysis using lysing solutions.
- Flow cytometer may produce false results if the device has not been aligned, calibrated and maintained appropriately.
- Data may be incorrectly interpreted if fluorescent signals were compensated wrongly or if gates were positioned inaccurately.
- Blood samples from some patients may exhibit abnormal values of positive cells.

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