## KOMBITEST™ CD3 FITC / CD19 PE

Cat. No. / Kat. č.: ED7053

### ENGLISH

#### 1. Specification

Content	50 tests, 1 ml	
Usage	20 µl per test (100 µl PBL)	
Specificity	CD3	CD19
Clone	UCHT1	4G7
Isotype (mouse)	IgG1	IgG1
Fluorochrome	FITC	PE
λ excitation	488 nm	488 nm
Emission maximum	525 nm	575 nm

#### 2. Intended use

The KOMBITEST™ CD3 FITC / CD19 PE is designed for percentage enumeration of mature human T and B lymphocytes in erythrocyte-lysed whole blood using Flow Cytometry.

#### 3. Principle

This test is based on the specific binding of monoclonal antibodies to the antigenic determinants expressed on the surface of leukocytes. The monoclonal antibodies are labeled with different fluorochromes which are excited via laser beam from a flow cytometer during analysis. Subsequent emissions of light from the fluorochromes of each cell are collected and analyzed by a flow cytometer. The fluorescence intensity differences enable the separation of cell subsets based on the expression of analyzed antigens.

The specific staining of blood cells is performed by the incubation of blood samples with the reagent followed by a lysis of red blood cells. Afterwards, unaffected leukocytes are subjected to analysis by a flow cytometer.

#### 4. Reagent provided

The reagent contains a premixed combination of mouse monoclonal antibody against human CD3 antigen (clone UCHT1) labeled with Fluorescein isothiocyanate (FITC), and mouse monoclonal antibody against human CD19 antigen (clone 4G7) labeled with R-phycoerythrin (PE). Labeled antibodies are diluted at optimum concentration in stabilizing phosphate buffered saline (PBS) solution containing 15mM sodium azide. The content of a vial (1 ml) is sufficient for 50 tests.

#### 5. Storage

Store vial in the dark at 2-8°C. Do not freeze.

#### 6. Precautions

- Intended for In Vitro Diagnostic use in laboratories outside USA and Canada. This CE-IVD reagent is in conformity with the European In Vitro Diagnostic Medical Device Directive 98/79/EC.
- Do not use after expiration date stamped on vial label.
- Avoid prolonged exposure to light.
- The content of the vial must not freeze.
- Avoid contamination of the reagent.
- The flow cytometer should be calibrated on a routine basis using fluorescent microbeads to ensure stable sensitivity of detectors.
- Any non-performance of the staining protocol may produce false results.
- The reagent contains sodium azide (NaN<sub>3</sub>) which is highly toxic in pure form. However, the concentration in the reagent (15mM) is not considered as hazardous. When disposing the reagent, flush the sink with a large volume of water.
- Blood samples are considered as potentially infectious and must be handled with care. Avoid all contact of the sample with the skin, eyes and mucosa.

#### 7. Necessary material not supplied

Test tubes for staining of blood samples (e.g. 12 × 75 mm), automatic pipettes with disposable tips, vortex mixer, commercial lysing solution, PBS buffer, flow cytometer.

#### 8. Blood specimen

The peripheral blood in a sterile tube with an anticoagulant (Heparin or EDTA) must be stored at room temperature and stained within 48 hours of drawing.

#### 9. Staining protocol

- Add 20 µl of KOMBITEST™ CD3 FITC / CD19 PE reagent to a test tube.
- Add 100 µl of blood sample to the tube. Vortex the tube.
- Incubate the tube for 15-20 minutes at room temperature in the dark.
- Perform lysis of red cells using lysing solution. It is recommended to use a commercial lysing solution containing formaldehyde as a fixative. Follow the instructions of the lysing solution manufacturer.
- Centrifuge tubes for 5 minutes at 300 g.
- Remove supernatant and resuspend pellet with 3-4 ml of PBS.
- Centrifuge tubes for 5 minutes at 300 g.
- Remove supernatant and resuspend pellet with 0.1 – 0.5 ml of PBS.

- Analyze the sample immediately using a flow cytometer or store sample at 2-8°C in the dark and analyze within 24 hours provided that cells were fixed.

## 10. Flow cytometry

Analyze stained samples using a flow cytometer equipped with excitation laser 488 nm and proper filters. Compensate fluorescent signals prior to or after data acquisition.

## 11. Data analysis

Visualize compensated data on the side-scatter (SSC) versus forward-scatter (FSC) plot. Set the gate for lymphocyte population as shown in figure 1. Alternatively, set the optimal lymphocyte gate using KOMBITEST™ CD45 FITC / CD14 PE (refer the datasheet for lymphocyte gate assessment procedure). Then make a CD3 FITC versus CD19 PE dot-plot of lymphocyte population as shown in figure 2. Separate populations using appropriate gate and calculate the percentage of T lymphocytes situated in lower-right quadrant (CD3+CD19-) and B lymphocytes situated in upper-left quadrant (CD3-CD19+) on the dot-plot.

## 12. Expected values

Results obtained in different laboratories may vary. Each laboratory should establish a normal range of cell subsets using its own test conditions. Results obtained in our laboratory are given in the table below.

Lymphocyte Subset	Unit	n	Mean	95% Range
CD3+CD19-	%	108	71	52-83
	cells/μl	54	1480	886-2341
CD3-CD19+	%	108	13	5-24
	cells/μl	54	242	87-445

## 13. Performance characteristics

### Specificity

The antibody UCHT1 recognizes the CD3 antigen of the TCR/CD3 complex on mature human T cells. The UCHT1 antibody reacts with the epsilon chain of the CD3 complex.


HLDA I; WS Code T 3  
 HLDA III; WS Code T 126  
 HLDA III; WS Code T 471  
 HLDA VI; WS Code T 6T-CD3.1

The mouse monoclonal antibody 4G7 recognizes CD19 (B4), a 95 kDa type I transmembrane glycoprotein of immunoglobulin superfamily, expressed on B lymphocytes and follicular dendritic cells; it is lost on plasma cells.


### Accuracy

The accuracy of the method was studied by the comparison of KOMBITEST™ with competitor's product in parallel staining of 43 blood samples. The regression analysis is given below.

### Regression Analysis of CD3+CD19- Lymphocytes


### Regression Analysis of CD3-CD19+ Lymphocytes


### Linearity

The linearity of the method was determined on 10 serial dilutions of leukocyte-enriched blood sample (buffy coat). Cell samples were stained by KOMBITEST™ in triplicates. Measured and expected values were expressed in terms of absolute count (cells/μl) in graphs given below.

### Range of CD3+CD19- Lymphocytes

### Range of CD3-CD19+ Lymphocytes

### Repeatability

The repeatability of the assay was measured on stabilized blood sample (Immuno-Troll™ Cells, Beckman-Coulter) in ten tubes in parallel. Coefficients of variation (CV) are given in the table below.

Lymphocyte Subset	Unit	n	Average	SD	CV
CD3+CD19-	%	10	70.3	0.94	1.34
CD3-CD19+	%	10	14.5	1.36	9.37

### Reproducibility

The reproducibility of the assay was measured on stabilized blood sample (Immuno-Troll™ Cells, Beckman-Coulter) under the same conditions for four weeks. Coefficients of variation (CV) are given in the table below.

Lymphocyte Subset	Unit	n	Average	SD	CV
CD3+CD19-	%	14	72.7	0.97	1.34
CD3-CD19+	%	14	13.2	1.16	8.78

### 14. Limitations

- Flow cytometer may produce false results if the device has not been aligned and maintained appropriately.
- Data may be incorrectly interpreted if fluorescent signals were compensated wrongly or if gates were positioned inaccurately.
- Blood samples from abnormal patients may exhibit abnormal values of positive cells.
- In case of a hyperleukocytose sample, it is recommended to dilute the blood sample with PBS to obtain leukocyte density approximately  $5 \times 10^6$  leukocytes/ml.

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
Fig. 1: Delimitation of lymphocyte population

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Barclay, Brown et al.: The Leukocyte Antigen FactsBook, 2nd edition, Harcourt Brace & Company, London, (1997).

Dubois B, Massacrier C, Caux C: Selective attraction of naive and memory B cells by dendritic cells. J Leukoc Biol. 2001 Oct;70(4):633-41.

Huang Y, Wange RL: T cell receptor signaling: beyond complex complexes. J Biol Chem. 2004 Jul 9;279(28):28827-30.

Kuhns MS, Davis MM, Garcia KC: Deconstructing the form and function of the TCR/CD3 complex. Immunity. 2006 Feb;24(2):133-9.

Leukocyte Typing III., McMichael A. J. et al (Eds.), Oxford University Press (1987).


Leukocyte Typing VI., Kishimoto T. et al. (Eds.), Garland Publishing Inc. (1997).

Rieux-Laucat F, Hivroz C, Lim A, Mateo V, Pellier I, Selz F, Fischer A, Le Deist F: Inherited and somatic CD3zeta mutations in a patient with T-cell deficiency. N Engl J Med. 2006 May 4;354(18):1913-21.

Siegers GM, Swamy M, Fernández-Malavé E, Minguet S, Rathmann S, Guardo AC, Pérez-Flores V, Regueiro JR, Alarcón B, Fisch P, Schamel WW: Different composition of the human and the mouse gammadelta T cell receptor explains different phenotypes of CD3gamma and CD3delta immunodeficiencies. J Exp Med. 2007 Oct 29;204(11):2537-44.

Torres PS, Alcover A, Zapata DA, Arnaud J, Pacheco A, Martín-Fernández JM, Villasevil EM, Sanal O, Regueiro JR: TCR dynamics in human mature T lymphocytes lacking CD3 gamma. J Immunol. 2003 Jun 15;170(12):5947-55.

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