T-cell BlastoFlowEx[®] Kit

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REF ED7642

ENGLISH

For Research Use Only. Not for use in diagnostic procedures.

EXBIO Praha, a.s.
Nad Safinou II 341
252 42 Vestec, Czech Republic
Tel: +420 261 090 666
Fax: +420 261 090 660
E-mail: orders@exbio.cz
www.exbio.cz

1. Recommended use

T-cell BlastoFlowEx[®] Kit is designed to measure the proliferative response of T-lymphocytes in activated samples of whole human blood. The kit uses anti-CD3 and anti-Ki67 antibody cocktail to detect the proliferating lymphocytes.

2. Components

T-cell BlastoFlowEx[®] Kit (ED7642) contains:

- ED7642-1: Fix and Lysing Solution, 10x concentrated, 1 x 25 ml,
- ED7642-2: Permeabilizing Solution, 10x concentrated, 1 x 10 ml,
- ED7642-3: EDTA solution, 1 x 2,5 ml,
- ED7642-4: CD3 APC/Ki-67 PE antibody cocktail, 2 x 2,5 ml.

Antibody cocktail specification

Content	2 x 2,5 ml	
Usage	50 µl per test	
Specificity	CD3	Ki-67
Clone	MEM-57	Ki-67
Isotype (mouse)	IgG2a	IgG1
Fluorochrome	APC	R-PE
λ excitation	633 nm	488 nm
Emission maximum	660 nm	575 nm

NOTE: The Fix and Lysing Solution and the Permeabilizing Solution must be diluted before use with deionized water (mix 1 volume of the concentrated solution with 9 volumes of deionized water).

The content of the kit is sufficient for 100 staining reactions.

3. Storage and shelf life after first opening

Reagents are stable until expiration date printed on the product label when stored at 2-8 °C.

Do not freeze the reagents.

Avoid prolonged exposure to light.

Shelf life after first opening is not different from shelf life printed on the product label.

Evidence of deterioration

If data obtained show any performance alteration, please contact manufacturer using the following e-mail address:

technical@exbio.cz

4. Warnings and precautions

Fix and Lysing Solution and **Permeabilizing Solution** contain 10% formaldehyde. The solutions are classified as hazardous according to the Regulation (EC) No 1272/2008. Wear protective gloves, protective clothing, eye protection and face protection when working with the reagents.

For detailed information about relevant H and P phrases please refer to the T-cell BlastoFlowEx[®] Kit Safety Data Sheet.

5. Principle of examination method

Anticoagulated blood samples are stimulated with mitogens for several days. At the end of incubation the sample contains a certain proportion of proliferating cells that show the so-called blastogenic transformation when viewed under a microscope. The kit uses measurements of intracellular Ki-67 protein expression, a proliferation specific marker, within CD3+ cells to detect the activated T-lymphocytes.

First the cells are incubated shortly in the presence of EDTA; the reagent stops intracellular signaling and proliferation by chelating calcium ions. Cells are then washed, fixed with formaldehyde, permeabilized with a non-ionic detergent containing buffer and stained with fluorescently labelled monoclonal antibodies. The proportion of Ki67 positive T-lymphocytes is analyzed with a flow cytometer.

6. Additional required materials and equipment

Materials required but not supplied

Deionized water

Phosphate buffered saline (PBS)

Automatic pipettes with disposable tips

Vortex mixer

Centrifuge with a rotor for 5ml test tubes (12 × 75 mm)

Waste container with disinfectant to collect the supernatants after cell centrifugations

Materials needed for blood stimulation

Cell culture medium X-VIVO™ (Lonza Cat. No. BE04-380Q)

Cell incubator (37 °C / 5% CO₂)

A laminar flow cabinet or laminar flow closet or tissue culture hood

Mitogens or antigens for stimulation

NOTE: The kit was optimized for use with whole blood cultures in 12 x 75 mm cytometry tubes. Single use cytometry tubes containing appropriate amounts of lyophilized mitogens are available from EXBIO under catalogue numbers:

ED7634 Phytohemagglutinin (PHA),

ED7635 Pokeweed mitogen (PWM),

ED7636 Concanavalin A (Con A),

ED7637 Stimulation Negative Control.

(please refer to EXBIO website for other consumables compatible with T-cell BlastoFlowEx® Kit)

Flow cytometer requirements

Flow cytometer equipped with a suitable excitation source and a set of detectors for light scatters and appropriate fluorescence (see antibody cocktail specification in section 2). The cytometer must be equipped with appropriate software for data acquisition and analysis. Laser alignment and regular cytometer calibration using fluorescent particles is highly advised.

7. Reagent preparation

- Specimen:

Blood samples should be collected into tubes containing **heparin anticoagulant** (Na-heparin or Li-heparin). Blood anticoagulated with EDTA or citrate is not suitable for stimulation.

It is recommended to use freshly withdrawn samples. Cell viability decreases in time after collection and may directly influence results.

WARNING: Human biological samples and blood specimens and any materials coming into contact with them are always considered as biohazard.

- Fix and Lysing Solution:

Prepare **2 ml** of the diluted solution for each staining reaction.

- Permeabilizing solution:

Prepare **0,5 ml** of the diluted solution for each staining reaction.

8. Procedure

Whole blood stimulation and cell culture

1. Place mitogen containing culture tubes inside a laminar flow cabinet.
(if using the EXBIO tubes, remove the caps, but keep them in the flow cabinet).
2. Add 500 µl of cell culture medium to each tube.
Add 50 µl of anticoagulated blood.
Close the tubes with caps and vortex gently to mix the content.
3. Ensure that all caps are in the ventilation position and place them in a cell incubator (37 °C with 5% CO₂) for at least 3 days.

Detection of Ki67 expressing T-lymphocytes

1. Remove the tubes from the incubator.
Remove and discard the caps.
Add **25 µl of EDTA** solution into each tube and mix.
2. Incubate for 10 minutes at 37 °C.
3. Add 2 ml of PBS, mix.
Centrifuge the cells for 5 minutes at 400g.
Decant supernatant.
4. Shake the tubes a little to disturb the pellet.
Add **2 ml of the diluted Fix and Lysing Solution**, mix.
Incubate for 10 minutes at room temperature.
5. Centrifuge the cells for 5 minutes at 400g.
Decant supernatant.
6. Add **0,5 ml of the diluted Permeabilizing Solution**, mix.
Incubate for 10 minutes at room temperature.
7. Add 2 ml of PBS.
Centrifuge the cells for 5 minutes at 400g.
Decant supernatant.
8. Add **50 µl of CD3/Ki-67** antibody cocktail, mix.
Incubate for at least 30 minutes at room temperature in the dark.
9. Add 2 ml of PBS.
Centrifuge the cells for 5 minutes at 400g.
Decant supernatant.
10. Resuspend the cells in 0,1-0,2 ml of 1% formaldehyde in PBS or Fix and Lysing solution diluted in PBS (mix 1 part of Fixation Buffer with 9 parts of PBS).
Store the processed samples at 2-8 °C in the dark until analysis.

9. Flow cytometry analysis

Please refer to "Instrument requirements" in section 6 for more information about flow cytometer requirements. Acquire data using software intended for sample analysis supplied with cytometer.

General cytometer setting recommendations:

- Set the threshold on forward light scatter so that only cells of interest are recorded and most of the debris excluded.
- Set the voltage on light scatter detectors, forward light scatter and side (perpendicular) light scatter so that the events of interest are on scale (consider that stimulated and unstimulated samples differ in their properties).
- Set voltage on fluorescence detectors so that all events of interest are on scale.
- Adjust the forward scatter area scaling factor (BD instruments) to enhance singlets vs doublets discrimination. When this factor is set properly the FSC peak area and the FSC peak height will have the same values and the singlet events will form a diagonal line (45 degrees, passing through zero).

Product specific notes:

- Acquire at least 3,000 of CD3+ lymphocytes per sample.
- Set the threshold on forward light scatter, do not set threshold on APC fluorescence detector. In our experience the threshold set on fluorescence caused distortion of results possibly due to the cytometer fluorescence signal/background processing algorithms.

Data analysis:

Look at the example analysis depicted at Figures 1 to 4 on page 8.

10. Trademarks

X-VIVO is a registered trademark of Lonza.

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Example data

Whole blood from a healthy blood donor was mixed 1:10 with X-VIVO™ culture medium and stimulated with phytohemagglutinin at concentration 20 µg/ml for 3 days and controls were incubated with X-VIVO™ medium only.

Plot the ungated events with side scatter on the Y-axis vs CD3 APC fluorescence on the X-axis. Draw a gate around **CD3+ events (Figure 1)**. Plot the events from CD3+ gate as the light forward scatter signal peak height vs light forward scatter signal peak area. Draw a diagonal gate around **singlets (Figure 2)**. Check the distribution of Ki-67 signal by displaying the CD3+/singlets as a dot plot with CD3 APC signal on the Y-axis vs Ki-67 PE signal on the X-axis. Look for the best **discrimination value** between the negative and positive population (dashed line) (**Figure 3**). Display the CD3+/singlets in histogram of PE fluorescence. Draw the **discrimination line** between the negative and positive peaks. Apply to all tubes and evaluate the **Ki-67+ percentage** in the sample (**Figure 4**).

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