

# Monoclonal Antibody to CD14, FITC conjugated (CD14 FITC)

Cat.No. ED7028

## 1. Intended purpose

The reagent CD14 FITC permits identification and enumeration of cell populations expressing human CD14 antigen in whole blood using flow cytometry.

## 2. Test principle

This test is based on specific binding of monoclonal antibody to the antigenic determinant expressed on the surface of leukocytes. The monoclonal antibody is labeled with fluorochrome which is excited via laser beam from a flow cytometer during analysis. Subsequent emission of light from fluorochromes of each cell is collected and analyzed by a flow cytometer. The fluorescence intensity differences enable separation of cell subsets based on expression of analyzed antigen. Specific staining of blood cells is performed by incubation of blood samples with the reagent followed by a lysis of red blood cells. Afterwards, unaffected leukocytes are subjected to analysis by a flow cytometer.

## 3. Reagents provided

The reagent contains mouse monoclonal antibody against human CD14 antigen (clone MEM-15) which was purified by affinity chromatography and labeled with Fluorescein isothiocyanate (FITC). The labeled antibody is diluted at optimum concentration in stabilizing phosphate buffered saline (PBS) solution containing 15mM sodium azide. The content of a vial (2 ml) is sufficient for 100 tests.

## Product specification

Content	100 tests, 2 ml
Usage	20 µl per test
Specificity	Human CD14
Clone	MEM-15
Isotype	Mouse IgG1
Fluorochrome	FITC
λ excitation	488 nm
Emission maximum	525 nm

## 4. Materials required but not provided

Test tubes for blood staining (e.g. 12 × 75 mm)  
Commercial lysing solution  
Phosphate buffered saline (PBS)  
Isotype control antibody (mouse IgG1 FITC)

## 5. Equipment required

Automatic pipettes with disposable tips  
Vortex mixer  
Centrifuge  
Flow cytometer with excitation laser 488 nm and proper filters

## 6. Storage and handling

Store the vial at 2-8 °C. Keep away from sunlight. Do not freeze. Do not aliquote. Expiration date is stated on a vial label and on outer packaging.

## 7. Warnings, precautions and limitations of use

- Intended for In Vitro Diagnostic use in laboratories outside USA and Canada. This CE-IVD kit is in conformity with the European Directive 98/79/EC.
- Do not use reagent after expiration date.
- Avoid reagents contamination.
- Avoid prolonged exposure to light.
- The content of the vial must not freeze.
- Any non-performance of staining protocol may produce false results.
- The reagent contains sodium azide (NaN<sub>3</sub>) which is highly toxic in pure form. However, the concentration in the reagent (15mM) is not considered as hazardous. When disposing the reagent, flush the sink with a large volume of water.
- Blood samples are considered as potentially infectious and must be handled with care. Avoid all contact of the sample with the skin, eyes and mucosa.
- In case of hyperleukocytose sample, it is recommended to dilute blood sample with PBS to obtain leukocyte density approximately 5 × 10<sup>6</sup> leukocytes/ml.
- Blood samples from abnormal patients may exhibit abnormal values of positive cells.
- Data may be incorrectly interpreted if fluorescent signals were compensated wrongly or if gates were positioned inaccurately.
- Flow cytometer may produce false results if

the device has not been aligned and maintained appropriately.

- Red blood cells from abnormal patients may be resistant to lysis using lysing solutions.
- Blood samples should be stained and analyzed within 24 hours from the blood collection.

## 8. Specimen

Use the peripheral human blood in a sterile tube with an anticoagulant (Heparin or EDTA). Blood must be stored at room temperature.

## 9. Procedure


### Staining protocol

- Add 20 µl of CD14 FITC reagent to a test tube, and the necessary amount of isotype control to a control tube.
- Add 100 µl of blood sample to each tube. Vortex the tubes.
- Incubate tubes for 20-30 minutes at room temperature in the dark.
- Perform lysis of red cells using lysing solution. It is recommended to use a commercial lysing solution containing formaldehyde as a fixative. Follow the instructions of the lysing solution manufacturer.
- Centrifuge tubes for 5 minutes at 300 g.
- Remove supernatant and resuspend pellet with 3-4 ml of PBS.
- Centrifuge tubes for 5 minutes at 300 g.
- Remove supernatant and resuspend pellet with 0.3 – 0.5 ml of PBS.
- Analyze samples immediately using flow cytometer or store samples at 2-8 °C in the dark and analyze within 24 hours provided that cells were fixed.

### Flow Cytometric Analysis

Analyze the sample stained with CD14 FITC using a flow cytometer. Visualize recorded data using appropriate plot such as side-scatter (SSC) versus FITC intensity as shown in figure 1. The brightest population (CD14+) belongs to monocytes. Set suitable gates for analysis. The region corresponding to the negative population should be set up using control cells which were stained by isotype control antibody.

Fig. 1: Leukocytes stained with CD14 FITC reagent



## 10. Analytical performance

### Specificity

The antibody MEM-15 reacts with CD14, a 53-55 kDa GPI (glycosylphosphatidylinositol) linked membrane glycoprotein expressed on monocytes, macrophages and weakly on granulocytes; it is expressed by most tissue macrophages. This antibody also reacts with soluble forms of CD14 found in serum and in the urine of some nephrotic patients. The monoclonal antibody MEM-15 was assigned to CD14 during the Human Leukocyte Differentiation Antigen workshop (HLDA3 WS Code: M 252).

### 11. Clinical performance

N/A

### 12. References

Asai Y et al. (2007) Soluble CD14 Discriminates Slight Structural Differences between Lipid As That Lead to Distinct Host Cell Activation. J Immunol. 179: 7674-83.

Lodrup Carlsen KC and Granum B (2007) Soluble CD14: role in atopic disease and recurrent infections, including otitis media. Curr Allergy Asthma Rep. 7: 436-43.

Fernández-Real JM et al. (2003) CD14 monocyte receptor, involved in the inflammatory cascade, and insulin sensitivity. J Clin Endocrinol Metab. 88: 1780-4.

Juan TS et al. (1995) Identification of a domain in soluble CD14 essential for lipopolysaccharide (LPS) signaling but not LPS binding. J Biol Chem. 270: 17237-42.

Bazil V et al. (1986) Biochemical characterization of a soluble form of the 53-kDa monocyte surface antigen. Eur J Immunol. 16: 1583-9.

Leukocyte Typing VI., Kishimoto T. et al. (Eds.), Garland Publishing Inc. (1997).

Leukocyte Typing V., Schlossman S. et al. (Eds.), Oxford University Press (1995).

Leukocyte Typing IV., Knapp W. et al. (Eds.),

Oxford University Press (1989).

Leukocyte Typing III., McMichael A. J. et al (Eds.), Oxford University Press (1987).

### 13. Manufacturer

EXBIO Praha, a.s.  
Nad Safinou II 341  
25250 Vestec  
Czech Republic

[info@exbio.cz](mailto:info@exbio.cz)

[technical@exbio.cz](mailto:technical@exbio.cz)

[orders@exbio.cz](mailto:orders@exbio.cz)

[www.exbio.cz](http://www.exbio.cz)

### 14. Trademarks

N/A

### 15. Revision History

- Version 1, ED7028\_IFU\_v1  
Initial Release
- Version 2, ED7028\_IFU\_v2  
The address was changed: "Nad Safinou II 341"
- Version 3, ED7028\_IFU\_v3  
Precautions section was changed. "Intended for professional use only." – removed. "Intended for In Vitro Diagnostic use in laboratories outside USA and Canada. This CE-IVD reagent is in conformity with the European In Vitro Diagnostic Medical Device Directive 98/79/EC." – added.
- Version 4, ED7028\_IFU\_v4  
Precautions section was changed – stabilizing added, solution – added and "0.2% (w/v) high-grade protease free Bovine Serum Albumin (BSA) as a stabilizing agent" – removed.
- Version 5, ED7028\_IFU\_v5  
The company logo changed. IFU layout changed. "Keep away from sunlight." – added. "Blood must be stored at room temperature.." - added. Postal code changed: "25250 Vestec"



## Monoclonal Antibody to CD14, FITC conjugated (CD14 FITC)

**100 tests** | Cat.No. ED7028



### Instructions for Use

Version: ED7028\_IFU\_v5\_EN

Date of Issue: 24-02-2020

EN

### Symbols



Catalogue number



Batch code



Use-by date



Temperature limits



Keep away from sunlight



In vitro diagnostic medical device



CE marking of conformity



Consult instructions for use



Manufacturer

The product is intended for In Vitro Diagnostic Use. In vivo diagnostic or therapeutic applications are strictly forbidden.

Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at [www.exbio.cz](http://www.exbio.cz). EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s..