



	Technical Data Sheet	
Product	Streptavidin FITC	
Cat. Number/Size	EXB0041	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Format	FITC	
Application	FC (QC tested)	
Description	Streptavidin FITC is intended for detecting biotinylated antigens and biotinylated antibodies by flow cytometry.	
Contents	Streptavidin FITC - 1 vial (0.5 ml) of Ready-to-use solution. No reconstitution is necessary. Reagent is provided in stabilizing phosphate buffered saline (PBS) solution containing 15mM sodium azide.	
Application details	Flow cytometry: Recomend use: 2.5 µg/ml. Extracellular and intracellular staining.	
Excitation laser	blue (488 nm)	
Preparation	Purified recombinant protein is conjugated with fluorescein isothiocyanate (FITC) under optimum conditions. Unconjugated protein and free fluorochrome are removed by size-exclusion chromatography.	
Formulation	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s.
EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz.
EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.