



	Technical Data Sheet	
Product	Annexin V APC	
Cat. Number/Size	EXB0028	100 tests
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Format	APC	
Application	FC (QC tested)	
Description	Annexin V APC is intended for detection of early apoptotic cells by flow cytometry.	
Contents	Annexin V APC - 1 vial (0.5 ml) of Ready-to-use solution. No reconstitution is necessary. Reagent is provided in stabilizing phosphate buffered saline (PBS) solution containing 15mM sodium azide. Sufficient for 100 tests.	
Excitation laser	red (633 nm)	
Preparation	Purified recombinant protein is conjugated with allophycocyanin (APC) under optimum conditions. Unconjugated protein and free fluorochrome are removed by size-exclusion chromatography.	
Formulation	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.	
Images and References	www.exbio.cz	

Unless indicated otherwise, all products are For Research Use Only and not for diagnostic use. In vivo diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s. © 2019 EXBIO Praha, a.s. All rights reserved.