



11-606-C025

Monoclonal Antibody to CD106 Purified Antibody (0.025 mg)

Clone:	STA
Isotype:	Mouse IgG1
Specificity:	The mouse monoclonal antibody STA recognizes CD106 antigen (VCAM-1), a 100-110 kDa type I membrane protein of the immunoglobulin superfamily, a crucial mediator of leukocyte adhesion, and a costimulation molecule. HLDA V; WS Code A013
Immunogen:	Human DS6 T cell line
Species Reactivity:	Human
Application:	Flow Cytometry Recommended dilution: 4-6 µg/ml Positive control: TNF-alpha activated HUVEC cells Immunoprecipitation Immunohistochemistry (frozen sections) Application note: acetone fixation ELISA Application note: capture mAb for soluble CD106
Purity:	> 95% (by SDS-PAGE)
Purification:	Purified from ascites by protein-A affinity chromatography.
Concentration:	1 mg/ml
Storage Buffer:	Phosphate buffered saline (PBS) with 15 mM sodium azide, approx. pH 7.4
Storage / Stability:	Store at 2-8°C. Do not use after expiration date stamped on vial label. For long-term storage aliquot and store at -20°C. Avoid freeze/thaw cycles.
Expiration:	See vial label
Lot Number:	See vial label
Background:	CD106 / VCAM-1 (vascular cell adhesion molecule-1) is an Ig-like cell surface adhesion molecule binding VLA-4 integrin. VCAM-1 is a potent T cell costimulatory molecule taking part in their positive selection and survival, as well as in adhesion, transendothelial migration and activation of peripheral T cells. VCAM-1 is also involved in endothelial cell-cell contacts. Whereas VCAM-1 normally mediates leukocyte extravasation to sites of tissue inflammation, tumour cells can use overexpressed VCAM-1 to escape T cell immunity. Soluble form of VCAM-1 (sVCAM-1) is an inflammatory marker and can be used also in prognosis of subsequent cardiovascular events following acute coronary syndromes.

For laboratory research only, not for drug, diagnostic or other use.



Antibodies

References:

*van Wetering S, van den Berk N, van Buul JD, Mul FP, Lommerse I, Mous R, ten Klooster JP, Zwaginga JJ, Hordijk PL: VCAM-1-mediated Rac signaling controls endothelial cell-cell contacts and leukocyte transmigration. *Am J Physiol Cell Physiol.* 2003 Aug;285(2):C343-52.

*Wu TC: The role of vascular cell adhesion molecule-1 in tumor immune evasion. *Cancer Res.* 2007 Jul 1;67(13):6003-6.

*Postadzhiyan AS, Tzontcheva AV, Kehayov I, Finkov B: Circulating soluble adhesion molecules ICAM-1 and VCAM-1 and their association with clinical outcome, troponin T and C-reactive protein in patients with acute coronary syndromes. *Clin Biochem.* 2007 Sep 19

*Paessens LC, Singh SK, Fernandes RJ, van Kooyk Y: Vascular cell adhesion molecule-1 (VCAM-1) and intercellular adhesion molecule-1 (ICAM-1) provide co-stimulation in positive selection along with survival of selected thymocytes. *Mol Immunol.* 2008 Jan;45(1):42-8.

*Leukocyte Typing V., Schlossman S. et al. (Eds.), Oxford University Press (1995).

*Leca G, Mansur SE, Bensussan A: Expression of VCAM-1 (CD106) by a subset of TCR gamma delta-bearing lymphocyte clones. Involvement of a metalloprotease in the specific hydrolytic release of the soluble isoform. *J Immunol.* 1995 Feb 1;154(3):1069-77.

*Yen YT, Liao F, Hsiao CH, Kao CL, Chen YC, Wu-Hsieh BA: Modeling the early events of severe acute respiratory syndrome coronavirus infection in vitro. *J Virol.* 2006 Mar;80(6):2684-93.

For laboratory research only, not for drug, diagnostic or other use.