

T cell workshop 2026

“T cells & friends - the rich tapestry of Immunology”



Program – Friday 4th September



9:00

Institute of Clinical Immunology and Allergology, First Faculty of Medicine (Studničkova 7, Albertov, Prague)
Organizational introduction & Infomation about the venues

9:10

Transfer to the venues

Faculty of Science - Viničná Microscopy Core Facility (VMCF)	The Institute of Hematology and Blood Transfusion	Institute of Clinical Immunology and Allergology
Viničná 7, Prague 2	U Nemocnice 2094/1, Prague 2	Studničkova 7, Prague 2
Josef Janda, Aleš Efenberk Practical demonstration of dry reagents in Flow Cytometry: DryFlowEx TBNK+ MonoDR 9-color (RUO) measured with CytoFLEX flow cytometer	Jan Musil, Antonín Ptáček BD FACS Discover S8 imaging spectral sorter: combining spectral setup with image analysis for cell sorting of human PBMC	Riccardo Pasculli (Accela) More than 20 Years of Microcapillary Flow Cytometry: Cytex Guava's innovation Ondřej Pelák (Waters) Beyond conventional data analysis in FlowJo

9:30	Block 1	Block 1	Block 1	Block 1
11:30				

11:35	Block 2	Block 2	Block 2	Block 2
13:35				

13:40

Lunch break

14:25	Block 3	Block 3	Block 3
16:25			

SPONSORS & PARTNERS

