



	Technical Data Sheet	
Product	Anti-Hu CD100 Purified	
Cat. Number/Size	11-772-C025	0.025 mg
	11-772-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	CD100	
Clone	133-1C6	
Format	Purified	
Reactivity	Human	
Application	FC (QC tested), IP, ELISA	
Application details	Flow cytometry: Recommended dilution: 1-4 µg/ml.	
Isotype	Mouse IgM	
Specificity	The mouse monoclonal antibody 133-1C6 recognizes an extracellular epitope of CD100, an approximately 150 kDa (when reduced) semaphorin family member expressed mainly on lymphocytes, NK cells, monocytes/macrophages and granulocytes, but also on some non-hematopoietic cells.	
Other names	COLL4, SEMA4D, semaphorin 4D	
Workshop	HLDA VI: WS Code N-L026	
Immunogen	PHA stimulated human PBL	
Entrez Gene ID	10507	
Gene name	SEMA4D	
NCBI Full Gene Name	semaphorin 4D	
UniProt ID	Q92854	
Concentration	1 mg/ml	
Preparation	Purified by sequential steps of physicochemical fractionation (differential precipitation and solid-phase chromatography methods).	
Formulation	Tris buffered saline (TBS), pH 8.0, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.