



	Technical Data Sheet	
Product	Anti-Hu CD57 Purified	
Cat. Number/Size	11-158-C025	0.025 mg
	11-158-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	CD57	
Clone	TB01	
Format	Purified	
Reactivity	Human	
Application	FC (QC tested), WB, IHC-P, IHC-F, ICC	
Application details	Flow cytometry: Recommended dilution: 1-4 µg/ml. Immunohistochemistry: Recommended dilution: 5-10 µg/ml.	
Isotype	Mouse IgM	
Specificity	The mouse monoclonal antibody TB01 recognizes CD57, a carbohydrate extracellular antigen present mainly on NK cells, NK T cells, and in neural tissue.	
Other names	HNK1, LEU7, 3-O-sulfoglucuronic acid beta1-3 Gal beta1-4 GlcNAc	
Workshop	HLDA VI: WS Code NK16	
Immunogen	A pool of neuroblastoma cell lines	
Entrez Gene ID	27087	
Gene name	B3GAT1	
NCBI Full Gene Name	beta-1,3-glucuronyltransferase 1	
UniProt ID	Q9P2W7	
Concentration	1 mg/ml	
Preparation	Purified by sequential steps of physicochemical fractionation (differential precipitation and solid-phase chromatography methods).	
Formulation	Tris buffered saline (TBS), pH 8.0, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.