



	Technical Data Sheet	
Product	Anti-Hu CD98 Purified	
Cat. Number/Size	11-268-C025	0.025 mg
	11-268-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	CD98	
Clone	MEM-108	
Format	Purified	
Reactivity	Human	
Application	FC (QC tested), IP, IHC-P	
Application details	Flow cytometry: Recommended dilution: 4 µg/ml. Immunohistochemistry (paraffin sections): Recommended dilution: 10 µg/ml; positive tissue: liver.	
Isotype	Mouse IgG1	
Specificity	The antibody MEM-108 reacts with an extracellular epitope of CD98, a 125 kDa disulfide-linked heterodimer (80 kDa glycosylated heavy chain + 45 kDa non-glycosylated light chain). CD98 is expressed on T lymphocytes (upon activation) and activated NK cells; it is also present at low levels on B lymphocytes, NK cells, monocytes and platelets.	
Other names	4F2, MDU1, 4F2HC, 4T2HC, NACAE, SLC3A2	
Workshop	HLDA VI: WS Code BP 409; HLDA VI: WS Code NL N-L017	
Immunogen	RAJI human Burkitt's lymphoma cell line	
Entrez Gene ID	8140	
Gene name	SLC7A5	
NCBI Full Gene Name	solute carrier family 7 member 5	
UniProt ID	Q01650	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.