

Technical Data Sheet

Product **Anti-Hu Ig Lambda Light Chain Purified**

Cat. Number/Size **11-317-C025** **0.025 mg**  
**11-317-C100** **0.1 mg**

**For Research Use Only.**

**Not for use in diagnostic or therapeutic procedures.**

---

Antigen	Ig Lambda Light Chain
Clone	Rs4
Format	Purified
Reactivity	Human
Application	ICC, ELISA
Isotype	Mouse IgG1
Specificity	The antibody Rs4 reacts with lambda light chains (22.5 kDa) of human immunoglobulin.
Other names	Immunoglobulin lambda, Igl

---

Concentration	1 mg/ml
Preparation	Purified by precipitation and chromatography
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze. Do not use after expiration date stamped on the label.

---

Images and References [www.exbio.cz](http://www.exbio.cz)

---

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at [www.exbio.cz](http://www.exbio.cz). EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.