



	Technical Data Sheet	
Product	Anti-NTAL Purified	
Cat. Number/Size	11-486-C025	0.025 mg
	11-486-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	NTAL	
Clone	PAb (486)	
Format	Purified	
Reactivity	Mouse	
Negative species	Human	
Application	WB	
Application details	Western blotting: Recommended dilution: 1 µg/ml; positive control: murine spleen; non-reducing conditions.	
Isotype	Rabbit polyclonal	
Specificity	The polyclonal antibody recognizes N-terminal cytoplasmic domain of mouse Non-T cell Activation Linker (NTAL), also known as LAB (linker of activated B cells), a 25 kDa transmembrane adaptor protein present in membrane microdomains (rafts) of B cells, NK cells and myeloid cells.	
Other names	LAB, LAT2	
Immunogen	Recombinant N-terminal cytoplasmic domain (aa 124-203) of mouse NTAL.	
Entrez Gene ID	7462	
Gene name	LAT2	
NCBI Full Gene Name	linker for activation of T cells family member 2	
UniProt ID	Q9GZY6	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.