

Technical Data Sheet

Product **Anti-Hu IgA Purified**  
Cat. Number/Size **11-947-C100** **0.1 mg**  
**For Research Use Only.**  
**Not for use in diagnostic or therapeutic procedures.**

---

Antigen	IgA
Clone	AD2.947
Format	Purified
Reactivity	Human
Application	IHC(P), ELISA
Isotype	Mouse IgG1
Specificity	The mouse monoclonal antibody AD2.947 reacts with human IgA1 and IgA2.
Other names	Immunoglobulin A
Immunogen	human IgA

---

Concentration	1 mg/ml
Preparation	Purified by sequential steps of physicochemical fractionation (differential precipitation and solid-phase chromatography methods).
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

---

Images and References [www.exbio.cz](http://www.exbio.cz)

---

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at [www.exbio.cz](http://www.exbio.cz). EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.