



	Technical Data Sheet	
Product	Anti-Hu deltaNp73 Purified	
Cat. Number/Size	11-889-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	deltaNp73	
Clone	DNp73-1.1	
Format	Purified	
Reactivity	Human	
Application	WB, ICC, IP	
Application details	Western blotting: Recommended dilution: 1 µg/ml	
Isotype	Mouse IgG1	
Specificity	The mouse monoclonal antibody DNp73-1.1 recognizes human deltaNp73; target epitope is MLYVGDPARHLAT.	
Other names	p73 DN	
Immunogen	KLH-coupled synthetic peptide MLYVGDPARHLATC	
Entrez Gene ID	7161	
Gene name	TP73	
NCBI Full Gene Name	tumor protein p73	
UniProt ID	O15350	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.