



|                       |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |                  |
|-----------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------|
|                       | Technical Data Sheet                                                                                                                                                                                                                                                                                                                                                                                                                                                                |                  |
| Product               | <b>Anti-Hu CD8 PE</b>                                                                                                                                                                                                                                                                                                                                                                                                                                                               |                  |
| Cat. Number/Size      | <b>1P-207-T025</b>                                                                                                                                                                                                                                                                                                                                                                                                                                                                  | <b>25 tests</b>  |
|                       | <b>1P-207-T100</b>                                                                                                                                                                                                                                                                                                                                                                                                                                                                  | <b>100 tests</b> |
|                       | <b>For Research Use Only.</b>                                                                                                                                                                                                                                                                                                                                                                                                                                                       |                  |
|                       | <b>Not for use in diagnostic or therapeutic procedures.</b>                                                                                                                                                                                                                                                                                                                                                                                                                         |                  |
| Antigen               | CD8                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |                  |
| Clone                 | MEM-31                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |                  |
| Format                | PE                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |                  |
| Reactivity            | Human                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |                  |
| Application           | FC (QC tested)                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |                  |
| Application details   | Flow cytometry: The reagent is designed for analysis of human blood cells using 20 µl reagent / 100 µl of whole blood or 10 <sup>6</sup> cells in a suspension. The content of a vial (2 ml) is sufficient for 100 tests.                                                                                                                                                                                                                                                           |                  |
| Excitation laser      | blue (488 nm)                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |                  |
| Isotype               | Mouse IgG2a                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |                  |
| Specificity           | The antibody MEM-31 recognizes a conformationally-dependent extracellular epitope of CD8, a cell surface glycoprotein found on most cytotoxic T lymphocytes that mediates efficient cell-cell interactions within the immune system. CD8 is a disulfide-linked dimer and exists as a CD8 alpha/alpha homodimer or CD8 alpha/beta heterodimer (each monomer approx. 32-34 kDa). The antibody does not react with formaldehyde-fixed cells; negative in Western blotting application. |                  |
| Other names           | p32, LEU2                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |                  |
| Workshop              | HLDA III: WS Code T 575                                                                                                                                                                                                                                                                                                                                                                                                                                                             |                  |
| Immunogen             | Crude thymus membrane fraction.                                                                                                                                                                                                                                                                                                                                                                                                                                                     |                  |
| Preparation           | Purified antibody is conjugated with R-phycoerythrin (PE) under optimum conditions. Unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.                                                                                                                                                                                                                                                                                                       |                  |
| Formulation           | Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide                                                                                                                                                                                                                                                                                                                                                                                                             |                  |
| Storage and handling  | Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.                                                                                                                                                                                                                                                                                                                                                                                                            |                  |
| Images and References | www.exbio.cz                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |                  |

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.