



	Technical Data Sheet	
Product	Anti-Hu C3aR Purified	
Cat. Number/Size	11-882-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	C3aR	
Clone	HC3aRZ8	
Format	Purified	
Reactivity	Human	
Application	FC (QC tested), IHC-F	
Application details	Flow cytometry: Recommended dilution: 1-4 µg/ml.	
Isotype	Mouse IgG2b	
Specificity	The mouse monoclonal antibody HC3aRZ8 recognizes an extracellular epitope of C3aR, a transmembrane chemotactic receptor for C3a anaphylatoxin.	
Other names	AZ3B, HNFAG09	
Immunogen	Human C3aR transfectants	
Entrez Gene ID	719	
Gene name	C3AR1	
NCBI Full Gene Name	complement C3a receptor 1	
UniProt ID	Q16581	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.