

Technical Data Sheet

Product

Anti-c-Myc tag Biotin

Cat. Number/Size

1B-433-C025

0.025 mg

1B-433-C100

0.1 mg

For Research Use Only.

Not for use in diagnostic or therapeutic procedures.

Antigen	c-Myc tag
Clone	9E10
Format	Biotin
Reactivity	Tagged fusion proteins in all species, Human
Application	FC (QC tested), IP, WB, IHC(P)
Application details	Flow cytometry: Intracellular or extracellular staining, depending on particular expression. Recommended dilution: 1-5 µg/ml.
Isotype	Mouse IgG1
Specificity	The antibody 9E10 can be used to detect the c-Myc tag.
Other names	bHLH, MRTL, MYCC
Immunogen	Synthetic peptide sequence (AEEQKLISEEDLL) corresponding to the C-terminal region of human c-Myc.
Entrez Gene ID	4609
Gene name	MYC
NCBI Full Gene Name	MYC proto-oncogene
UniProt ID	P01106

Concentration	1 mg/ml
Preparation	Purified antibody is conjugated with biotin LC-NHS ester under optimum conditions and unconjugated antibody and free biotin are removed by size-exclusion chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.