

Technical Data Sheet

Product	Anti-LRRC32 Purified	
Cat. Number/Size	11-744-C025	0.025 mg
	11-744-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	

Antigen	LRRC32
Clone	GARP5
Format	Purified
Reactivity	Human
Application	WB
Application details	Western blotting: Recommended dilution: 1-2 µg/ml; positive control: human thrombocytes, reducing conditions
Isotype	Mouse IgG1
Specificity	The mouse monoclonal antibody GARP5 recognizes GARP / LRRC32, an approximately 80 kDa glycoprotein expressed e.g. on the surface of megakaryocytes, platelets and activated Treg cells.
Other names	GARP, D11S833E
Immunogen	Purified human sGARP protein
Entrez Gene ID	2615
Gene name	LRRC32
NCBI Full Gene Name	leucine rich repeat containing 32
UniProt ID	Q14392

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.