

Technical Data Sheet

Product	Anti-Cytokeratins Biotin	
Cat. Number/Size	1B-108-C025	0.025 mg
	1B-108-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	

Antigen	Cytokeratins
Clone	C-11
Format	Biotin
Reactivity	Mammalian
Application	WB, FC (QC tested), IP, IHC(P), ICC
Isotype	Mouse IgG1
Specificity	The antibody C-11 reacts with Cytokeratin peptides 4, 5, 6, 8, 10, 13, 18. Cytokeratins are a member of intermediate filaments subfamily represented in epithelial tissues.
Other names	cytokeratin, CYK, CK, KRT
Immunogen	Keratin-enriched preparation from human epidermoid carcinoma cell line A431.

Concentration	1 mg/ml
Preparation	The purified antibody is conjugated with biotin-LC-NHS under optimum conditions. The reagent is free of unconjugated biotin.
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze. Do not use after expiration date stamped on vial label.

Images and References	www.exbio.cz
-----------------------	--

Unless indicated otherwise, all products are For Research Use Only and not for diagnostic use. In vivo diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s. © 2019 EXBIO Praha, a.s. All rights reserved.