

Technical Data Sheet

Product	Anti-PilT (F. tularensis ssp. tularensis) Purified	
Cat. Number/Size	11-180-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	

Antigen	PilT (F. tularensis ssp. tularensis)
Clone	PAb (180)
Format	Purified
Reactivity	Francisella tularensis ssp. tularensis
Application	WB, ELISA
Isotype	Rabbit polyclonal
Specificity	The rabbit polyclonal antibody to PilT reacts with Francisella tularensis subsp. tularensis antigen PilT (FTT0088, also known as FTT_0088). The antibody does not react with Francisella tularensis subsp. holarctica homolog.
Other names	Type IV pili NBP
Immunogen	Recombinant PilT protein produced in E. coli
Entrez Gene ID	3190704
Gene name	pilT
NCBI Full Gene Name	Type IV pili nucleotide-binding protein
UniProt ID	Q5NII0

Concentration	1 mg/ml
Preparation	Purified from rabbit serum by affinity chromatography
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze. Do not use after expiration date stamped on the label.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.