

Technical Data Sheet

Product	<b>Anti-Hu C5aR2 Purified</b>	
Cat. Number/Size	<b>11-880-C100</b>	<b>0.1 mg</b>
	<b>For Research Use Only.</b>	
	<b>Not for use in diagnostic or therapeutic procedures.</b>	

---

Antigen	C5aR2
Clone	1D9-M12
Format	Purified
Reactivity	Human, Non-human primates
Application	WB, FC (QC tested)
Application details	Flow cytometry: Recommended dilution: 2-6 µg/ml
Isotype	Mouse IgG2a kappa
Specificity	The mouse monoclonal antibody 1D9-M12 recognizes an extracellular epitope on C5aR2 (C5L2), a C5a complement receptor, which is coexpressed with C5aR1 (CD88) in neutrophils, as well as e.g. in mast cells, astrocytes, or macrophages.
Other names	C5L2, GPR77, GPF77
Immunogen	L1.2 cells transfected with human C5aR2
Entrez Gene ID	27202
Gene name	C5AR2
NCBI Full Gene Name	complement component 5a receptor 2
UniProt ID	Q9P296

---

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

---

Images and References [www.exbio.cz](http://www.exbio.cz)

---

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at [www.exbio.cz](http://www.exbio.cz). EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.