



	Technical Data Sheet	
Product	Anti-Hu CD62P Biotin	
Cat. Number/Size	1B-667-C025	0.025 mg
	1B-667-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	CD62P	
Clone	AK4	
Format	Biotin	
Reactivity	Human, Non-human primates	
Application	FC (QC tested), IP, WB, ICC, IHC-F	
Application details	Flow cytometry: Recommended dilution: 2-5 µg/ml.	
Isotype	Mouse IgG1 kappa	
Specificity	The antibody AK4 recognizes an extracellular epitope of CD62P (P-selectin), a 140 kD single chain type I transmembrane glycoprotein present in secretory alpha-granules in platelets, in Weibel-Palade bodies in endothelial cells and in megakaryocytes; it is relocated to the plasma membrane upon activation.	
Other names	P-selectin, GMP140, SELP, PADGEM, GRMP, LECAM3, PSEL	
Workshop	HLDA VI: WS Code P-44	
Immunogen	Human platelets	
Entrez Gene ID	6403	
Gene name	SELP	
NCBI Full Gene Name	selectin P	
UniProt ID	P16109	
Concentration	1 mg/ml	
Preparation	Purified antibody is conjugated with biotin LC-NHS ester under optimum conditions and unconjugated antibody and free biotin are removed by size-exclusion chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.