



	Technical Data Sheet	
Product	Anti-Hu CD194 Purified	
Cat. Number/Size	11-961-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	CD194	
Clone	rHIQB002	
Format	Purified	
Reactivity	Human	
Application	FC (QC tested)	
Application details	Flow cytometry: Recommended dilution: 0.5-4 µg/ml.	
Isotype	Mouse IgG1	
Specificity	The recombinant mouse monoclonal antibody rHIQB002 recognizes an extracellular epitope of human chemokine receptor CD194 (CCR4).	
Other names	CCR4, CC-CKR-4, CKR4, CMKBR4, ChemR13, HGCN	
Entrez Gene ID	1233	
Gene name	CCR4	
NCBI Full Gene Name	C-C motif chemokine receptor 4	
UniProt ID	P51679	
Concentration	1 mg/ml	
Preparation	Purified by affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.