

Technical Data Sheet

Product	Anti-Hu CD27 Purified	
Cat. Number/Size	11-308-C025	0.025 mg
	11-308-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	

Antigen	CD27
Clone	LT27
Format	Purified
Reactivity	Human
Application	FC (QC tested)
Application details	Flow cytometry: Recommended dilution: 2 µg/ml.
Isotype	Mouse IgG2a
Specificity	The antibody LT27 reacts with an extracellular epitope of CD27 (T14), a 50-55 kDa type I transmembrane glycoprotein (member of the TNF-receptor superfamily) expressed on medullary thymocytes, peripheral T lymphocytes, some B lymphocytes and NK cells.
Other names	T14, TNFRSF7, LPFS2, S152, Tp55
Workshop	HLDA V: WS Code T T-CD27.01
Immunogen	Human peripheral blood lymphocytes
Entrez Gene ID	939
Gene name	CD27
NCBI Full Gene Name	CD27 molecule
UniProt ID	P26842

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.