

Technical Data Sheet

Product	Anti-Hu CD226 Purified	
Cat. Number/Size	11-926-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	

Antigen	CD226
Clone	11A8
Format	Purified
Reactivity	Non-human primates, Human
Application	FC (QC tested)
Application details	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Isotype	Mouse IgG1 kappa
Specificity	The mouse monoclonal antibody 11A8 recognizes an extracellular epitope within CD226 transmembrane glycoprotein.
Other names	DNAM1, PTA1, TLISA1
Entrez Gene ID	10666
Gene name	CD226
NCBI Full Gene Name	CD226 molecule
UniProt ID	Q15762

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.