



	Technical Data Sheet	
Product	Anti-Hu SUSD2 APC	
Cat. Number/Size	1A-835-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	SUSD2	
Clone	W5C5	
Format	APC	
Reactivity	Human	
Application	FC (QC tested)	
Application details	Flow cytometry: Recommended dilution: 1-4 µg/ml	
Excitation laser	red (633 nm)	
Isotype	Mouse IgG1	
Specificity	The mouse monoclonal antibody W5C5 recognizes an extracellular epitope of SUSD2, a type I transmembrane protein expressed on mesenchymal stem-like cells. This antibody selectively binds to a MSCs in both bone marrow and endometrium or tonsil, and can be used for their identification and isolation.	
Immunogen	WERI-RB-1 retinoblastoma cells	
Entrez Gene ID	56241	
Gene name	SUSD2	
NCBI Full Gene Name	sushi domain containing 2	
UniProt ID	Q9UGT4	
Concentration	0.1 mg/ml	
Preparation	Purified antibody is conjugated with activated allophycocyanin (APC) under optimum conditions and unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.	
Formulation	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.