

Technical Data Sheet

Product	<b>Anti-Hu CD335 Purified</b>	
Cat. Number/Size	<b>11-853-C100</b>	<b>0.1 mg</b>
	<b>For Research Use Only.</b>	
	<b>Not for use in diagnostic or therapeutic procedures.</b>	

---

Antigen	CD335
Clone	9E2
Format	Purified
Reactivity	Human
Application	FC (QC tested)
Application details	Flow cytometry: Recommended dilution: 1-5 µg/ml
Isotype	Mouse IgG1
Specificity	The mouse monoclonal antibody 9E2 recognizes an extracellular epitope of human CD335, a transmembrane glycoprotein, serving as a natural cytotoxicity receptor against virus-infected or cancer cells.
Other names	NKP46, NCR1, LY94
Immunogen	fusion protein of human CD335 and Fc part of IgG
Entrez Gene ID	9437
Gene name	NCR1
NCBI Full Gene Name	natural cytotoxicity triggering receptor 1
UniProt ID	O76036

---

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

---

Images and References [www.exbio.cz](http://www.exbio.cz)

---

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at [www.exbio.cz](http://www.exbio.cz). EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.