



	Technical Data Sheet	
Product	Anti-Hu GARP Purified	
Cat. Number/Size	11-744-C025	0.025 mg
	11-744-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	GARP	
Clone	GARP5	
Format	Purified	
Reactivity	Human	
Application	WB, FC	
Application details	Western blotting: Recommended dilution: 1-2 µg/ml; positive control: human thrombocytes, reducing conditions	
Isotype	Mouse IgG1	
Specificity	The mouse monoclonal antibody GARP5 recognizes GARP / LRRC32, an approximately 80 kDa glycoprotein expressed e.g. on the surface of megakaryocytes, platelets and activated Treg cells.	
Other names	LRRC32, garpin, D11S833E	
Immunogen	Purified human sGARP protein	
Entrez Gene ID	2615	
Gene name	LRRC32	
NCBI Full Gene Name	leucine rich repeat containing 32	
UniProt ID	Q14392	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.