

Technical Data Sheet

Product

**Anti-Hu CD19 FITC**

Cat. Number/Size

**1F-663-T025**

**25 tests**

**1F-663-T100**

**100 tests**

**For Research Use Only.**

**Not for use in diagnostic or therapeutic procedures.**

Antigen	CD19
Clone	4G7
Format	FITC
Reactivity	Human
Application	FC (QC tested)
Application details	Flow cytometry: The reagent is designed for analysis of human blood cells using 20 µl reagent / 100 µl of whole blood or 10 <sup>6</sup> cells in a suspension. The content of a vial (2 ml) is sufficient for 100 tests.
Excitation laser	blue (488 nm)
Isotype	Mouse IgG1 kappa
Specificity	The mouse monoclonal antibody 4G7 recognizes an extracellular epitope of human CD19.
Other names	B4, Leu-12, CVID3
Workshop	HLDA: WS Code 2 B43
Immunogen	Human CCL (chronic lymphocytic leukemia) cells
Entrez Gene ID	930
Gene name	CD19
NCBI Full Gene Name	CD19 molecule
UniProt ID	P15391

Preparation	Purified antibody is conjugated with fluorescein isothiocyanate (FITC) under optimum conditions and unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.
Formulation	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.

Images and References [www.exbio.cz](http://www.exbio.cz)

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at [www.exbio.cz](http://www.exbio.cz). EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.