



|                       |                                                                                                                                    |               |
|-----------------------|------------------------------------------------------------------------------------------------------------------------------------|---------------|
|                       | Technical Data Sheet                                                                                                               |               |
| Product               | <b>Anti-Hu CD85h Purified</b>                                                                                                      |               |
| Cat. Number/Size      | <b>11-925-C100</b>                                                                                                                 | <b>0.1 mg</b> |
|                       | <b>For Research Use Only.</b>                                                                                                      |               |
|                       | <b>Not for use in diagnostic or therapeutic procedures.</b>                                                                        |               |
| Antigen               | CD85h                                                                                                                              |               |
| Clone                 | 24                                                                                                                                 |               |
| Format                | Purified                                                                                                                           |               |
| Reactivity            | Non-human primates, Human                                                                                                          |               |
| Application           | FC (QC tested), IHC-F                                                                                                              |               |
| Application details   | Flow cytometry: Recommended dilution: 1-4 µg/ml.                                                                                   |               |
| Isotype               | Mouse IgG2b kappa                                                                                                                  |               |
| Specificity           | The mouse monoclonal antibody 24 recognizes an extracellular epitope of CD85h, an approximately 70 kDa transmembrane glycoprotein. |               |
| Other names           | ILT1, LIR7, LILRA2                                                                                                                 |               |
| Immunogen             | human CD85h transfectants                                                                                                          |               |
| Entrez Gene ID        | 11027                                                                                                                              |               |
| Gene name             | LILRA2                                                                                                                             |               |
| NCBI Full Gene Name   | leukocyte immunoglobulin like receptor A2                                                                                          |               |
| UniProt ID            | Q8N149                                                                                                                             |               |
| Concentration         | 1 mg/ml                                                                                                                            |               |
| Preparation           | Purified by protein-A affinity chromatography.                                                                                     |               |
| Formulation           | Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide                                                                        |               |
| Storage and handling  | Store at 2-8°C. Do not freeze.                                                                                                     |               |
| Images and References | <a href="http://www.exbio.cz">www.exbio.cz</a>                                                                                     |               |

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at [www.exbio.cz](http://www.exbio.cz). EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.