



	Technical Data Sheet	
Product	Anti-DDIT4L Purified	
Cat. Number/Size	11-725-C025	0.025 mg
	11-725-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	DDIT4L	
Clone	DDIT-03	
Format	Purified	
Reactivity	Human	
Application	WB, FC-IC	
Application details	Flow cytometry: Recommended dilution: 1-4 µg/ml. Intracellular staining.	
Isotype	Mouse IgG1	
Specificity	The mouse monoclonal antibody DDIT-03 recognizes DDIT4L / REDD2 protein, which belongs to intracellular stress-induced proteins involved in mediation of cell death.	
Other names	DDIT4L, REDD2, Rtp801L	
Immunogen	N-terminal recombinant fragment of human DDIT4L (amino acids 2-98)	
Entrez Gene ID	115265	
Gene name	DDIT4L	
NCBI Full Gene Name	DNA damage inducible transcript 4 like	
UniProt ID	Q96D03	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.