

Technical Data Sheet

Product	Anti-PDE8A Purified	
Cat. Number/Size	11-710-C025	0.025 mg
	11-710-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	

Antigen	PDE8A
Clone	EM-52
Format	Purified
Reactivity	Human
Application	WB
Application details	Western blotting: Recommended dilution: 1-2 µg/ml.
Isotype	Mouse IgG1 kappa
Specificity	The mouse monoclonal antibody EM-52 recognizes PDE8a (phosphodiesterase 8 A), an approximately 93 kDa intracellular enzyme that hydrolyzes cAMP to 5'AMP.
Other names	HsT19550
Immunogen	Human PDE8a
Entrez Gene ID	5151
Gene name	PDE8A
NCBI Full Gene Name	phosphodiesterase 8A
UniProt ID	O60658

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.