

Technical Data Sheet

Product **Anti-Hu CD138 Purified Azide Free**

Cat. Number/Size **10-520-C025** **0.025 mg**  
**10-520-C100** **0.1 mg**

**For Research Use Only.**

**Not for use in diagnostic or therapeutic procedures.**

---

Antigen	CD138
Clone	B-A38
Format	Purified Azide Free
Reactivity	Human
Application	FC (QC tested), IHC
Isotype	Mouse IgG1
Specificity	The antibody B-A38 recognizes CD138 (syndecan 1), a 65-70 kDa heparan sulfate proteoglycan expressed mainly in the epidermis and plasma cells, but also in growth factor-stimulated lymphocytes.
Other names	syndecan 1, SYND1, SDC
Immunogen	U266 human peripheral blood myeloma cell line

---

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography
Formulation	Azide free phosphate buffered saline (PBS), approx. pH 7.4; 0.2 µm filter sterilized
Storage and handling	Store at 2-8°C. Do not freeze. Do not use after expiration date stamped on vial label.

---

Images and References [www.exbio.cz](http://www.exbio.cz)

---

Unless indicated otherwise, all products are For Research Use Only and not for diagnostic use. In vivo diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at [www.exbio.cz](http://www.exbio.cz). EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s. © 2019 EXBIO Praha, a.s. All rights reserved.