



	Technical Data Sheet	
Product	Anti-Hu CD169 Purified	
Cat. Number/Size	11-797-C025	0.025 mg
	11-797-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	CD169	
Clone	7-239	
Format	Purified	
Reactivity	Human	
Application	FC (QC tested), IP, WB, IHC-F	
Application details	Flow cytometry: Recommended dilution: 1-4 µg/ml	
Isotype	Mouse IgG1	
Specificity	The mouse monoclonal antibody 7-239 recognizes an extracellular epitope of CD169 (sialoadhesin, Siglec-1), a 210 kDa type I transmembrane glycoprotein expressed on macrophages and dendritic cells.	
Other names	SIGLEC1, SN, sialoadhesin	
Immunogen	human rhinovirus 14-infected monocyte-derived dendritic cells	
Entrez Gene ID	6614	
Gene name	SIGLEC1	
NCBI Full Gene Name	sialic acid binding Ig like lectin 1	
UniProt ID	Q9BZZ2	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.