

Technical Data Sheet

Product	Anti-OPAL1 Purified	
Cat. Number/Size	11-699-C025	0.025 mg
	11-699-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	

Antigen	OPAL1
Clone	OPAL1-01
Format	Purified
Reactivity	Human
Application	FC (QC tested), WB
Application details	Flow cytometry: Recommended dilution: 1-4 µg/ml, intracellular staining. Western blotting: Recommended dilution: 1-2 µg/ml, non-reducing conditions.
Isotype	Mouse IgG2a
Specificity	The mouse monoclonal antibody OPAL1-01 recognizes an intracellular epitope of OPAL1 (outcome predictor in acute leukemia 1), a transmembrane adaptor protein expressed mainly by megacaryocytes and platelets.
Other names	WBP1L
Immunogen	Recombinant fragment of human OPAL1 (amino acids 152-342)
Entrez Gene ID	54838
Gene name	WBP1L
NCBI Full Gene Name	WW domain binding protein 1 like
UniProt ID	Q9NX94

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Tris buffered saline (TBS), pH 8.0, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.