

Technical Data Sheet

Product	<b>Anti-DR3 Purified</b>	
Cat. Number/Size	<b>11-757-C025</b>	<b>0.025 mg</b>
	<b>11-757-C100</b>	<b>0.1 mg</b>
	<b>For Research Use Only.</b>	
	<b>Not for use in diagnostic or therapeutic procedures.</b>	

---

Antigen	DR3
Clone	JD3
Format	Purified
Reactivity	Human
Application	FC (QC tested)
Isotype	Mouse IgG1
Specificity	The mouse monoclonal antibody JD3 recognizes an extracellular epitope of DR3 (APO-3, TNFRSF25), a transmembrane protein of TNFR superfamily expressed mainly in lymphocyte-enriched tissues.
Other names	TNFRSF25, TR3, DDR3, LARD, APO-3, TRAMP, WSL-1, WSL-LR, TNFRSF12
Immunogen	human DR3-Ig fusion protein
Entrez Gene ID	8718
Gene name	TNFRSF25
NCBI Full Gene Name	TNF receptor superfamily member 25
UniProt ID	Q93038

---

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze. Do not use after expiration date stamped on the label.

---

Images and References [www.exbio.cz](http://www.exbio.cz)

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at [www.exbio.cz](http://www.exbio.cz). EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.