

Technical Data Sheet

Product	Anti-HLA-G Purified	
Cat. Number/Size	11-436-C025	0.025 mg
	11-436-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	

Antigen	HLA-G
Clone	MEM-G/2
Format	Purified
Reactivity	Human
Application	IHC-P, WB, IHC-F
Application details	Western blotting: Recommended dilution 1-2 µg/ml.
Isotype	Mouse IgG1
Specificity	The antibody MEM-G/2 recognizes an extracellular epitope on free heavy chain of all the HLA-G isoforms. HLA-G belongs to the MHC Class I molecules (MHC Class Ib; nonclassical) and it is expressed on the surface of trophoblast cells.
Immunogen	Bacterially expressed recombinant HLA-G heavy chain (denatured).
Entrez Gene ID	3135
Gene name	HLA-G
NCBI Full Gene Name	major histocompatibility complex, class I, G
UniProt ID	P17693

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.