



	Technical Data Sheet	
Product	Anti-Ms CD11a FITC	
Cat. Number/Size	1F-622-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	CD11a	
Clone	M17/4	
Format	FITC	
Reactivity	Mouse	
Application	FC (QC tested)	
Application details	Flow cytometry: Recommended dilution: 1-2 µg/ml.	
Excitation laser	blue (488 nm)	
Isotype	Rat IgG2a	
Specificity	The rat monoclonal antibody M17/4 reacts with an extracellular epitope of CD11a (alpha-subunit of murine LFA-1), a 180 kDa type I transmembrane glycoprotein expressed on B and T lymphocytes, monocytes, macrophages, neutrophils, basophils and eosinophils.	
Other names	LFA-1, LFA1A, ITGAL	
Immunogen	C57BL/6 mouse splenic secondary cytotoxic T lymphocytes	
Entrez Gene ID	16408	
Gene name	Itgal	
NCBI Full Gene Name	integrin subunit alpha L	
UniProt ID	E9Q5M7	
Concentration	0.5 mg/ml	
Preparation	Purified antibody is conjugated with fluorescein isothiocyanate (FITC) under optimum conditions and unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.