



Technical Data Sheet

Product	Anti-Hu IgE Purified	
Cat. Number/Size	11-324-C025	0.025 mg
	11-324-C100	0.1 mg
For Research Use Only.		
Not for use in diagnostic or therapeutic procedures.		
Antigen	IgE	
Clone	BE5	
Format	Purified	
Reactivity	Human	
Application	FC (QC tested), ELISA	
Application details	Flow cytometry: Recommended dilution: 5 µg/ml.	
Isotype	Mouse IgG1	
Specificity	The antibody BE5 reacts with human IgE; it recognizes an epitope different from the ones recognized by 4G7.325 and 4H10 antibodies to IgE.	
Other names	Immunoglobulin E	
Immunogen	Purified human IgE.	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.