

Technical Data Sheet

Product	<b>Anti-Ms CD26 Purified</b>	
Cat. Number/Size	<b>11-776-C025</b>	<b>0.025 mg</b>
	<b>11-776-C100</b>	<b>0.1 mg</b>
	<b>For Research Use Only.</b>	
	<b>Not for use in diagnostic or therapeutic procedures.</b>	

---

Antigen	CD26
Clone	H194-112
Format	Purified
Reactivity	Mouse
Application	IP, FC (QC tested), IHC-F
Application details	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Isotype	Rat IgG2a kappa
Specificity	The rat monoclonal antibody H194-112 recognizes an extracellular epitope of CD26, a 110 kDa type II membrane glycoprotein, which is a peptidase expressed on mature thymocytes, T cells (especially activated), B cells, NK cells and macrophages.
Other names	DPP4, TP103, ADCP-2, ADABP
Immunogen	BALB/c thymocytes
Entrez Gene ID	13482
Gene name	Dpp4
NCBI Full Gene Name	dipeptidyl peptidase 4
UniProt ID	P28843

---

Concentration	1 mg/ml
Preparation	Purified by protein-G affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

---

Images and References [www.exbio.cz](http://www.exbio.cz)

---

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at [www.exbio.cz](http://www.exbio.cz). EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.